

Package ‘MSstatsTMT’

May 11, 2024

Title Protein Significance Analysis in shotgun mass spectrometry-based proteomic experiments with tandem mass tag (TMT) labeling

Version 2.12.0

Date 2024-04-23

Description The package provides statistical tools for detecting differentially abundant proteins in shotgun mass spectrometry-based proteomic experiments with tandem mass tag (TMT) labeling. It provides multiple functionalities, including aata visualization, protein quantification and normalization, and statistical modeling and inference. Furthermore, it is inter-operable with other data processing tools, such as Proteome Discoverer, MaxQuant, OpenMS and SpectroMine.

License Artistic-2.0

Depends R (>= 4.2)

Imports limma, lme4, lmerTest, methods, data.table, stats, utils, ggplot2, grDevices, graphics, MSstats, MSstatsConvert, checkmate

Suggests BiocStyle, knitr, rmarkdown, testthat

VignetteBuilder knitr

biocViews ImmunoOncology, MassSpectrometry, Proteomics, Software

Encoding UTF-8

LazyData true

URL <http://msstats.org/msstatstmt/>

BugReports <https://groups.google.com/forum/#!forum/msstats>

RoxygenNote 7.2.3

git_url <https://git.bioconductor.org/packages/MSstatsTMT>

git_branch RELEASE_3_19

git_last_commit 200b85c

git_last_commit_date 2024-04-30

Repository Bioconductor 3.19

Date/Publication 2024-05-10

Author Devon Kohler [aut, cre],
 Ting Huang [aut],
 Meena Choi [aut],
 Mateusz Staniak [aut],
 Sicheng Hao [aut],
 Olga Vitek [aut]

Maintainer Devon Kohler <kohler.d@northeastern.edu>

Contents

| | |
|---|----|
| .checkContrastMatrix | 3 |
| .checkSummarizationParams | 3 |
| .countRunsWithNorm | 5 |
| .documentFunction | 5 |
| .getNormalizationAbundance | 6 |
| .getPhilosopherInput | 7 |
| .getRunsMedian | 7 |
| .handleSingleContrastTMT | 8 |
| .logSum | 8 |
| .logSummarizationParams | 9 |
| .makeContrastSingleTMT | 10 |
| .makeFactorColumnsTMT | 10 |
| .medianPolish | 11 |
| .normalizePeptides | 11 |
| .normalizeProteins | 12 |
| .prepareForSummarization | 12 |
| .removeRedundantChannels | 13 |
| .summarizeMSstats | 13 |
| .summarizeSimpleStat | 14 |
| .summarizeTMP | 15 |
| .summarizeTMT | 15 |
| annotation.mine | 16 |
| annotation.mq | 17 |
| annotation.pd | 18 |
| dataProcessPlotsTMT | 19 |
| evidence | 21 |
| getProcessedTMT | 22 |
| getSummarizedTMT | 22 |
| groupComparisonTMT | 23 |
| input.pd | 24 |
| MaxQtoMSstatsTMTFormat | 25 |
| MSstatsComparisonModelSingleTMT | 27 |
| MSstatsFitComparisonModelsTMT | 27 |
| MSstatsGroupComparisonOutputTMT | 28 |
| MSstatsGroupComparisonTMT | 28 |
| MSstatsModerateTTest | 29 |
| MSstatsNormalizeTMT | 29 |

MSstatsPrepareForGroupComparisonTMT 30

MSstatsPrepareForSummarizationTMT 30

MSstatsSummarizationOutputTMT 31

MSstatsSummarizeTMT 32

MSstatsTestSingleProteinTMT 33

MSstatsTMT 33

OpenMStoMSstatsTMTFormat 34

PDtoMSstatsTMTFormat 35

PhilosophertoMSstatsTMTFormat 37

proteinGroups 38

proteinSummarization 39

quant.pd.msstats 41

raw.mine 42

raw.om 42

raw.pd 43

SpectroMinetoMSstatsTMTFormat 44

test.pairwise 46

Index **47**

.checkContrastMatrix *check whether pairwise comparison. If pairwise, generate a contrast matrix.*

Description

check whether pairwise comparison. If pairwise, generate a contrast matrix.

Usage

.checkContrastMatrix(contrast_matrix)

Value

a contrast matrix

.checkSummarizationParams
Check validity of parameters to proteinSummarization function

Description

Check validity of parameters to proteinSummarization function

Usage

```
.checkSummarizationParams(
  data,
  method,
  global_norm,
  reference_norm,
  remove_norm_channel,
  remove_empty_channel,
  MBimpute,
  maxQuantileforCensored
)
```

Arguments

| | |
|-------------------------------------|---|
| <code>data</code> | Name of the output of PDtoMSstatsTMTFormat function or peptide-level quantified data from other tools. It should have columns ProteinName, PeptideSequence, Charge, PSM, Mixture, TechRepMixture, Run, Channel, Condition, BioReplicate, Intensity |
| <code>method</code> | Four different summarization methods to protein-level can be performed : "msstats"(default), "MedianPolish", "Median", "LogSum". |
| <code>global_norm</code> | Global median normalization on peptide level data (equalizing the medians across all the channels and MS runs). Default is TRUE. It will be performed before protein-level summarization. |
| <code>reference_norm</code> | Reference channel based normalization between MS runs on protein level data. TRUE(default) needs at least one reference channel in each MS run, annotated by 'Norm' in Condition column. It will be performed after protein-level summarization. FALSE will not perform this normalization step. If data only has one run, then reference_norm=FALSE. |
| <code>remove_norm_channel</code> | TRUE(default) removes 'Norm' channels from protein level data. |
| <code>remove_empty_channel</code> | TRUE(default) removes 'Empty' channels from protein level data. |
| <code>MBimpute</code> | only for method="msstats". TRUE (default) imputes missing values by Accelerated failure model. FALSE uses minimum value to impute the missing value for each peptide precursor ion. |
| <code>maxQuantileforCensored</code> | We assume missing values are censored. maxQuantileforCensored is Maximum quantile for deciding censored missing value, for instance, 0.999. Default is Null. |

Value

TRUE invisibly if all parameters are valid

.countRunsWithNorm *Utility function: count runs with "Norm" channel*

Description

Utility function: count runs with "Norm" channel

Usage

```
.countRunsWithNorm(run, condition)
```

Arguments

| | |
|-----------|----------------------------|
| run | vector of run labels |
| condition | vector of condition labels |

Value

integer

.documentFunction *A dummy function to store shared documentation items.*

Description

A dummy function to store shared documentation items.

Usage

```
.documentFunction(  
  fewMeasurements,  
  useUniquePeptide,  
  summaryforMultipleRows,  
  removeProtein_with1Feature,  
  removeProtein_with1Protein,  
  removeOxidationMpeptides,  
  removeMpeptides  
)
```

Arguments

| | |
|---|---|
| <code>fewMeasurements</code> | <code>'remove'</code> (default) will remove the features that have 1 or 2 measurements across runs. |
| <code>useUniquePeptide</code> | TRUE (default) removes peptides that are assigned for more than one proteins. We assume to use unique peptide for each protein. |
| <code>summaryForMultipleRows</code> | <code>max</code> (default) or <code>sum</code> - when there are multiple measurements for certain feature and certain run, use highest or sum of multiple intensities. |
| <code>removeProtein_with1Feature</code> | TRUE will remove the proteins which have only 1 feature, which is the combination of peptide, precursor charge, fragment and charge. FALSE is default. |
| <code>removeOxidationMpeptides</code> | TRUE will remove the peptides including 'oxidation (M)' in modification. FALSE is default. |
| <code>removeMpeptides</code> | TRUE will remove the peptides including 'M' sequence. FALSE is default. |
| <code>removeProtein_with1Peptide</code> | TRUE will remove the proteins which have only 1 peptide and charge. FALSE is default. |
| <code>use_log_file</code> | logical. If TRUE, information about data processing will be saved to a file. |
| <code>append</code> | logical. If TRUE, information about data processing will be added to an existing log file. |
| <code>verbose</code> | logical. If TRUE, information about data processing will be printed to the console. |
| <code>log_file_path</code> | character. Path to a file to which information about data processing will be saved. If not provided, such a file will be created automatically. If <code>'append = TRUE'</code> , has to be a valid path to a file. |

Value

NULL.

.getNormalizationAbundance

Utility function: get mean abundance for "Norm" channels

Description

Utility function: get mean abundance for "Norm" channels

Usage

`.getNormalizationAbundance(abundance, condition)`

Arguments

| | |
|-----------|----------------------------|
| abundance | vector of abundances |
| condition | vector of condition labels |

Value

numeric

`.getPhilosopherInput` *Convert Philosopher parameters to consistent format*

Description

Convert Philosopher parameters to consistent format

Usage

```
.getPhilosopherInput(input, path, folder)
```

Arguments

| | |
|-------|--|
| input | data.frame of 'msstats.csv' file produced by Philosopher |
|-------|--|

`.getRunsMedian` *Utility function: get median from unique values per run*

Description

Utility function: get median from unique values per run

Usage

```
.getRunsMedian(input)
```

Arguments

| | |
|-------|-------------------|
| input | data.table / list |
|-------|-------------------|

Value

numeric

```
.handleSingleContrastTMT
```

perform statistical inference for single protein and single contrast

Description

perform statistical inference for single protein and single contrast

Usage

```
.handleSingleContrastTMT(  
  contrast,  
  fit,  
  single_protein,  
  coefs,  
  protein,  
  groups,  
  s2_posterior,  
  rho,  
  vss,  
  df_prior,  
  s2_df  
)
```

```
.logSum
```

Utility function: compute log of sum of 2^x

Description

Utility function: compute log of sum of 2^x

Usage

```
.logSum(x)
```

Arguments

x numeric

Value

numeric

.logSummarizationParams

Log parameters for proteinSummarization function

Description

Log parameters for proteinSummarization function

Usage

```
.logSummarizationParams(  
  method,  
  global_norm,  
  reference_norm,  
  remove_norm_channel,  
  remove_empty_channel  
)
```

Arguments

| | |
|----------------------|---|
| method | Four different summarization methods to protein-level can be performed : "msstats"(default), "MedianPolish", "Median", "LogSum". |
| global_norm | Global median normalization on peptide level data (equalizing the medians across all the channels and MS runs). Default is TRUE. It will be performed before protein-level summarization. |
| reference_norm | Reference channel based normalization between MS runs on protein level data. TRUE(default) needs at least one reference channel in each MS run, annotated by 'Norm' in Condition column. It will be performed after protein-level summarization. FALSE will not perform this normalization step. If data only has one run, then reference_norm=FALSE. |
| remove_norm_channel | TRUE(default) removes 'Norm' channels from protein level data. |
| remove_empty_channel | TRUE(default) removes 'Empty' channels from protein level data. |

Value

TRUE invisibly after logging successfully

`.makeContrastSingleTMT`*Make a contrast*

Description

Make a contrast

Usage`.makeContrastSingleTMT(fit, contrast, single_protein, coefs)`**Value**

a contrast vector

`.makeFactorColumnsTMT` *Converts required columns to factor in summarization output*

Description

Converts required columns to factor in summarization output

Usage`.makeFactorColumnsTMT(input)`**Arguments**`input` `data.table`**Value**

a data table with factored columns

`.medianPolish` *Tukey median polish*

Description

Tukey median polish

Usage

```
.medianPolish(intensities, num_channels)
```

Arguments

`intensities` vector of log-intensities per protein and run
`num_channels` number of channels

Value

numeric vector with length 'num_channels'

`.normalizePeptides` *Normalization between channels (before summarization)*

Description

Normalization between channels (before summarization)

Usage

```
.normalizePeptides(input, normalize)
```

Arguments

`input` data.table
`normalize` logical, if TRUE, 'input' data will be normalized

Value

data.table

`.normalizeProteins` *Normalization between MS runs (after protein summarization)*

Description

Normalization between MS runs (after protein summarization)

Usage

```
.normalizeProteins(input, normalize)
```

Arguments

| | |
|-----------|---|
| input | data.table |
| normalize | logical, if TRUE, data will be normalized |

Value

data.table

`.prepareForSummarization`
Prepare TMT data for protein-level summarization

Description

Prepare TMT data for protein-level summarization

Usage

```
.prepareForSummarization(input)
```

Arguments

| | |
|-------|------------|
| input | data.table |
|-------|------------|

Value

data.table with required column types

`.removeRedundantChannels`
Remove empty and normalization channels

Description

Remove empty and normalization channels

Usage

```
.removeRedundantChannels(input, remove_empty_channel, remove_norm_channel)
```

Arguments

`input` data.table processed by the protein summarization function
`remove_empty_channel`
 TRUE(default) removes 'Empty' channels from protein level data.
`remove_norm_channel`
 TRUE(default) removes 'Norm' channels from protein level data.

Value

data.table

`.summarizeMSstats` *Summarization based on MSstats*

Description

Summarization based on MSstats

Usage

```
.summarizeMSstats(  
  input,  
  annotation,  
  impute,  
  max_quantile_censored = NULL,  
  log_file_path = NULL  
)
```

Arguments

| | |
|-----------------------|---|
| input | data.table |
| annotation | data.table with run and channel annotation |
| impute | only for method="msstats". TRUE (default) imputes missing values by Accelerated failure model. FALSE uses minimum value to impute the missing value for each peptide precursor ion. |
| max_quantile_censored | We assume missing values are censored. maxQuantileforCensored is Maximum quantile for deciding censored missing value, for instance, 0.999. Default is Null. |
| log_file_path | path to a MSstats log file |

Value

data.table

`.summarizeSimpleStat` *Summarize TMT data with a simple aggregate of log-intensities*

Description

Summarize TMT data with a simple aggregate of log-intensities

Usage

```
.summarizeSimpleStat(input, annotation, stat_aggregate)
```

Arguments

| | |
|----------------|---|
| input | data.table |
| annotation | data.table with run and channel annotation |
| stat_aggregate | function that will be used to compute protein-level summary |

Value

data.table

.summarizeTMP *Summarize TMT data with median polish*

Description

Summarize TMT data with median polish

Usage

```
.summarizeTMP(input, annotation)
```

Arguments

| | |
|------------|--|
| input | data.table |
| annotation | data.table with run and channel annotation |

Value

data.table with summarized protein intensities

.summarizeTMT *Performs summarization for TMT data*

Description

Performs summarization for TMT data

Usage

```
.summarizeTMT(  
  input,  
  method,  
  annotation,  
  impute,  
  max_quantile_censored,  
  log_file_path  
)
```

Arguments

| | |
|------------|--|
| input | data.table |
| method | "mstats"/"MedianPolish"/"LogSum"/"Median" |
| annotation | data.table with run and channel annotation |

| | |
|-----------------------|---|
| impute | only for method="msstats". TRUE (default) imputes missing values by Accelerated failure model. FALSE uses minimum value to impute the missing value for each peptide precursor ion. |
| max_quantile_censored | We assume missing values are censored. maxQuantileforCensored is Maximum quantile for deciding censored missing value, for instance, 0.999. Default is Null. |
| log_file_path | path to a MSstats log file |

Value

data.table

| | |
|-----------------|---|
| annotation.mine | <i>Example of annotation file for raw.mine, which is the output of SpectroMine.</i> |
|-----------------|---|

Description

Annotation of example data, raw.mine, in this package. It should be prepared by users. The variables are as follows:

Usage

annotation.mine

Format

A data frame with 72 rows and 7 variables.

Details

- Run : MS run ID. It should be the same as R.FileName info in raw.mine
- Channel : Labeling information (TMT6_126, ..., TMT6_131). The channels should be consistent with the channel columns in raw.mine.
- Condition : Condition (ex. Healthy, Cancer, Time0). If the channel doesn't have sample, please add 'Empty' under Condition.
- Mixture : Mixture of samples labeled with different TMT reagents, which can be analyzed in a single mass spectrometry experiment.
- TechRepMixture : Technical replicate of one mixture. One mixture may have multiple technical replicates. For example, if 'TechRepMixture' = 1, 2 are the two technical replicates of one mixture, then they should match with same 'Mixture' value.
- Fraction : Fraction ID. One technical replicate of one mixture may be fractionated into multiple fractions to increase the analytical depth. Then one technical replicate of one mixture should correspond to multiple fractions. For example, if 'Fraction' = 1, 2, 3 are three fractions of the first technical replicate of one TMT mixture of biological subjects, then they should have same 'TechRepMixture' and 'Mixture' value.

- BioReplicate : Unique ID for biological subject. If the channel doesn't have sample, please add 'Empty' under BioReplicate

Examples

```
head(annotation.mine)
```

| | |
|---------------|--|
| annotation.mq | <i>Example of annotation file for evidence, which is the output of MaxQuant.</i> |
|---------------|--|

Description

Annotation of example data, evidence, in this package. It should be prepared by users. The variables are as follows:

Usage

```
annotation.mq
```

Format

A data frame with 150 rows and 7 variables.

Details

- Run : MS run ID. It should be the same as Raw.file info in raw.mq
- Channel : Labeling information (channel.0, ..., channel.9). The channel index should be consistent with the channel columns in raw.mq.
- Condition : Condition (ex. Healthy, Cancer, Time0)
- Mixture : Mixture of samples labeled with different TMT reagents, which can be analyzed in a single mass spectrometry experiment. If the channel doesn't have sample, please add 'Empty' under Condition.
- TechRepMixture : Technical replicate of one mixture. One mixture may have multiple technical replicates. For example, if 'TechRepMixture' = 1, 2 are the two technical replicates of one mixture, then they should match with same 'Mixture' value.
- Fraction : Fraction ID. One technical replicate of one mixture may be fractionated into multiple fractions to increase the analytical depth. Then one technical replicate of one mixture should correspond to multiple fractions. For example, if 'Fraction' = 1, 2, 3 are three fractions of the first technical replicate of one TMT mixture of biological subjects, then they should have same 'TechRepMixture' and 'Mixture' value.
- BioReplicate : Unique ID for biological subject. If the channel doesn't have sample, please add 'Empty' under BioReplicate.

Examples

```
head(annotation.mq)
```

| | |
|---------------|--|
| annotation.pd | <i>Example of annotation file for raw.pd, which is the PSM output of Proteome Discoverer</i> |
|---------------|--|

Description

Annotation of example data, raw.pd, in this package. It should be prepared by users. The variables are as follows:

Usage

```
annotation.pd
```

Format

A data frame with 150 rows and 7 variables.

Details

- Run : MS run ID. It should be the same as Spectrum.File info in raw.pd.
- Channel : Labeling information (126, ... 131). It should be consistent with the channel columns in raw.pd.
- Condition : Condition (ex. Healthy, Cancer, Time0)
- Mixture : Mixture of samples labeled with different TMT reagents, which can be analyzed in a single mass spectrometry experiment. If the channel doesn't have sample, please add 'Empty' under Condition.
- TechRepMixture : Technical replicate of one mixture. One mixture may have multiple technical replicates. For example, if 'TechRepMixture' = 1, 2 are the two technical replicates of one mixture, then they should match with same 'Mixture' value.
- Fraction : Fraction ID. One technical replicate of one mixture may be fractionated into multiple fractions to increase the analytical depth. Then one technical replicate of one mixture should correspond to multiple fractions. For example, if 'Fraction' = 1, 2, 3 are three fractions of the first technical replicate of one TMT mixture of biological subjects, then they should have same 'TechRepMixture' and 'Mixture' value.
- BioReplicate : Unique ID for biological subject. If the channel doesn't have sample, please add 'Empty' under BioReplicate.

Examples

```
head(annotation.pd)
```

dataProcessPlotsTMT *Visualization for explanatory data analysis - TMT experiment*

Description

To illustrate the quantitative data and quality control of MS runs, dataProcessPlotsTMT takes the quantitative data and summarized data from function 'proteinSummarization' as input and generate two types of figures in pdf files as output : (1) profile plot (specify "ProfilePlot" in option type), to identify the potential sources of variation for each protein; (2) quality control plot (specify "QCPlot" in option type), to evaluate the systematic bias between MS runs and channels.

Usage

```
dataProcessPlotsTMT(
  data,
  type,
  ylimUp = FALSE,
  ylimDown = FALSE,
  x.axis.size = 10,
  y.axis.size = 10,
  text.size = 4,
  text.angle = 90,
  legend.size = 7,
  dot.size.profile = 2,
  ncol.guide = 5,
  width = 10,
  height = 10,
  which.Protein = "all",
  originalPlot = TRUE,
  summaryPlot = TRUE,
  address = ""
)
```

Arguments

| | |
|-------------|--|
| data | the output of proteinSummarization function. It is a list with data frames 'FeatureLevelData' and 'ProteinLevelData' |
| type | choice of visualization. "ProfilePlot" represents profile plot of log intensities across MS runs. "QCPlot" represents box plots of log intensities across channels and MS runs. |
| ylimUp | upper limit for y-axis in the log scale. FALSE(Default) for Profile Plot and QC Plot uses the upper limit as rounded off maximum of log ₂ (intensities) after normalization + 3.. |
| ylimDown | lower limit for y-axis in the log scale. FALSE(Default) for Profile Plot and QC Plot uses 0.. |
| x.axis.size | size of x-axis labeling for "Run" and "channel in Profile Plot and QC Plot. |

| | |
|-------------------------------|---|
| <code>y.axis.size</code> | size of y-axis labels. Default is 10. |
| <code>text.size</code> | size of labels represented each condition at the top of Profile plot and QC plot. Default is 4. |
| <code>text.angle</code> | angle of labels represented each condition at the top of Profile plot and QC plot. Default is 0. |
| <code>legend.size</code> | size of legend above Profile plot. Default is 7. |
| <code>dot.size.profile</code> | size of dots in Profile plot. Default is 2. |
| <code>ncol.guide</code> | number of columns for legends at the top of plot. Default is 5. |
| <code>width</code> | width of the saved pdf file. Default is 10. |
| <code>height</code> | height of the saved pdf file. Default is 10. |
| <code>which.Protein</code> | Protein list to draw plots. List can be names of Proteins or order numbers of Proteins. Default is "all", which generates all plots for each protein. For QC plot, "allonly" will generate one QC plot with all proteins. |
| <code>originalPlot</code> | TRUE(default) draws original profile plots, without normalization. |
| <code>summaryPlot</code> | TRUE(default) draws profile plots with protein summarization for each channel and MS run. |
| <code>address</code> | the name of folder that will store the results. Default folder is the current working directory. The other assigned folder has to be existed under the current working directory. An output pdf file is automatically created with the default name of "ProfilePlot.pdf" or "QCplot.pdf". The command address can help to specify where to store the file as well as how to modify the beginning of the file name. If address=FALSE, plot will be not saved as pdf file but showed in window. |

Value

plot or pdf

Examples

```
data(input.pd)
quant.msstats = proteinSummarization(input.pd,
                                     method="msstats",
                                     global_norm=TRUE,
                                     reference_norm=TRUE)

## Profile plot
dataProcessPlotsTMT(data=quant.msstats,
                    type='ProfilePlot',
                    width = 21,
                    height = 7)

## NottoRun: QC plot
# dataProcessPlotsTMT(data=quant.msstats,
#                     # type='QCPlot',
#                     # width = 21,
```

```
# height = 7)
```

evidence

Example of output from MaxQuant for TMT-10plex experiments.

Description

Example of evidence.txt from MaxQuant. It is the input for MaxQtoMSstatsTMTFormat function, with proteinGroups.txt and annotation file. Annotation file should be made by users. It includes peak intensities for 10 proteins among 15 MS runs with TMT10. The important variables are as follows:

Usage

```
evidence
```

Format

A data frame with 1075 rows and 105 variables.

Details

- Proteins
- Protein.group.IDs
- Modified.sequence
- Charge
- Raw.file
- Score
- Potential.contaminant
- Reverse
- Channels : Reporter.intensity.corrected.0, ..., Reporter.intensity.corrected.9

Examples

```
head(evidence)
```

getProcessedTMT *Get processed feature-level data*

Description

Get processed feature-level data

Usage

```
getProcessedTMT(summarized, input)
```

Arguments

summarized output of the MSstatsSummarizeTMT function
input output of MSstatsNormalizeTMT function

Value

data.table

getSummarizedTMT *Get protein-level data from MSstatsSummarizeTMT output*

Description

Get protein-level data from MSstatsSummarizeTMT output

Usage

```
getSummarizedTMT(summarized)
```

Arguments

summarized output of the MSstatsSummarizeTMT function

Value

data.table

| | |
|--------------------|---|
| groupComparisonTMT | <i>Finding differentially abundant proteins across conditions in TMT experiment</i> |
|--------------------|---|

Description

Tests for significant changes in protein abundance across conditions based on a family of linear mixed-effects models in TMT experiment. Experimental design of case-control study (patients are not repeatedly measured) is automatically determined based on proper statistical model.

Usage

```
groupComparisonTMT(
  data,
  contrast.matrix = "pairwise",
  moderated = FALSE,
  adj.method = "BH",
  remove_norm_channel = TRUE,
  remove_empty_channel = TRUE,
  save_fitted_models = FALSE,
  use_log_file = TRUE,
  append = FALSE,
  verbose = TRUE,
  log_file_path = NULL
)
```

Arguments

| | |
|----------------------|--|
| data | the output of proteinSummarization function. It is a list with data frames 'FeatureLevelData' and 'ProteinLevelData' |
| contrast.matrix | Comparison between conditions of interests. 1) default is "pairwise", which compare all possible pairs between two conditions. 2) Otherwise, users can specify the comparisons of interest. Based on the levels of conditions, specify 1 or -1 to the conditions of interests and 0 otherwise. The levels of conditions are sorted alphabetically. |
| moderated | TRUE will moderate t statistic; FALSE (default) uses ordinary t statistic. |
| adj.method | adjusted method for multiple comparison. "BH" is default. |
| remove_norm_channel | TRUE(default) removes "Norm" channels from protein level data. |
| remove_empty_channel | TRUE(default) removes "Empty" channels from protein level data. |
| save_fitted_models | logical, if TRUE, fitted models will be added to |
| use_log_file | logical. If TRUE, information about data processing will be saved to a file. |

| | |
|---------------|---|
| append | logical. If TRUE, information about data processing will be added to an existing log file. |
| verbose | logical. If TRUE, information about data processing will be printed to the console. |
| log_file_path | character. Path to a file to which information about data processing will be saved. If not provided, such a file will be created automatically. If 'append = TRUE', has to be a valid path to a file. |

Value

a list that consists of the following elements: (1) ComparisonResult: statistical testing results; (2) FittedModel: the fitted linear models

Examples

```
data(input.pd)
# use protein.summarization() to get protein abundance data
quant.pd.msstats = proteinSummarization(input.pd,
                                       method="msstats",
                                       global_norm=TRUE,
                                       reference_norm=TRUE)

test.pairwise = groupComparisonTMT(quant.pd.msstats, moderated = TRUE)
head(test.pairwise$ComparisonResult)

# Only compare condition 0.125 and 1
levels(quant.pd.msstats$ProteinLevelData$Condition)

# Compare condition 1 and 0.125
comparison=matrix(c(-1,0,0,1),nrow=1)

# Set the names of each row
row.names(comparison)="1-0.125"

# Set the column names
colnames(comparison)= c("0.125", "0.5", "0.667", "1")
test.contrast = groupComparisonTMT(data = quant.pd.msstats,
contrast.matrix = comparison,
moderated = TRUE)
head(test.contrast$ComparisonResult)
```

input.pd

Example of output from PDtoMSstatsTMTFormat function

Description

It is made from [raw.pd](#) and [annotation.pd](#), which is the output of PDtoMSstatsTMTFormat function. It should include the required columns as below.

Usage

```
input.pd
```

Format

A data frame with 20110 rows and 11 variables.

Details

- ProteinName : Protein ID
- PeptideSequence : peptide sequence
- Charge : peptide charge
- PSM : peptide ion and spectra match
- Channel : Labeling information (126, ... 131)
- Condition : Condition (ex. Healthy, Cancer, Time0)
- BioReplicate : Unique ID for biological subject.
- Run : MS run ID
- Mixture : Unique ID for TMT mixture.
- TechRepMixture : Unique ID for technical replicate of one TMT mixture.
- Intensity: Protein Abundance

Examples

```
head(input.pd)
```

MaxQtoMSstatsTMTFormat

Generate MSstatsTMT required input format from MaxQuant output

Description

Generate MSstatsTMT required input format from MaxQuant output

Usage

```
MaxQtoMSstatsTMTFormat(  
  evidence,  
  proteinGroups,  
  annotation,  
  which.proteinid = "Proteins",  
  rmProt_Only.identified.by.site = FALSE,  
  useUniquePeptide = TRUE,  
  rmPSM_withfewMea_withinRun = TRUE,
```

```

    rmProtein_with1Feature = FALSE,
    summaryforMultipleRows = sum,
    use_log_file = TRUE,
    append = FALSE,
    verbose = TRUE,
    log_file_path = NULL,
    ...
)

```

Arguments

| | |
|--------------------------------|--|
| evidence | name of 'evidence.txt' data, which includes feature-level data. |
| proteinGroups | name of 'proteinGroups.txt' data. |
| annotation | data frame which contains column Run, Fraction, TechRepMixture, Mixture, Channel, BioReplicate, Condition. Refer to the example 'annotation.mq' for the meaning of each column. |
| which.proteinid | Use 'Proteins' (default) column for protein name. 'Leading.proteins' or 'Leading.razor.proteins' or 'Gene.names' can be used instead to get the protein ID with single protein. However, those can potentially have the shared peptides. |
| rmProt_Only.identified.by.site | TRUE will remove proteins with '+' in 'Only.identified.by.site' column from proteinGroups.txt, which was identified only by a modification site. FALSE is the default. |
| useUniquePeptide | TRUE(default) removes peptides that are assigned for more than one proteins. We assume to use unique peptide for each protein. |
| rmPSM_withfewMea_withinRun | TRUE (default) will remove the features that have 1 or 2 measurements within each Run. |
| rmProtein_with1Feature | TRUE will remove the proteins which have only 1 peptide and charge. Default is FALSE. |
| summaryforMultipleRows | sum(default) or max - when there are multiple measurements for certain feature in certain run, select the feature with the largest summation or maximal value. |
| use_log_file | logical. If TRUE, information about data processing will be saved to a file. |
| append | logical. If TRUE, information about data processing will be added to an existing log file. |
| verbose | logical. If TRUE, information about data processing will be printed to the console. |
| log_file_path | character. Path to a file to which information about data processing will be saved. If not provided, such a file will be created automatically. If 'append = TRUE', has to be a valid path to a file. |
| ... | additional parameters to 'data.table::fread'. |

Value

data.frame of class "MSstatsTMT"

Examples

```
head(evidence)
head(proteinGroups)
head(annotation.mq)
input.mq <- MaxQtoMSstatsTMTFormat(evidence, proteinGroups, annotation.mq)
head(input.mq)
```

MSstatsComparisonModelSingleTMT

Fit a linear model for group comparison for a single protein

Description

Fit a linear model for group comparison for a single protein

Usage

```
MSstatsComparisonModelSingleTMT(single_protein, protein_name)
```

Arguments

`single_protein` protein-level data for a single protein (single element of list created by the `MSstatsPrepareForGroupComparisonTMT` function)

`protein_name` name of a protein from the `single_protein` data.table

Value

list

MSstatsFitComparisonModelsTMT

Fit linear models for group comparison

Description

Fit linear models for group comparison

Usage

```
MSstatsFitComparisonModelsTMT(input)
```

Arguments

input output of the MSstatsPrepareForGroupComparisonTMT function

Value

list

MSstatsGroupComparisonOutputTMT

Combine testing results for individual proteins

Description

Combine testing results for individual proteins

Usage

```
MSstatsGroupComparisonOutputTMT(testing_results, adj_method)
```

Arguments

testing_results

output of the MSstatsGroupComparisonTMT function

adj_method

method that will be used to adjust p-values for multiple comparisons

Value

data.table

MSstatsGroupComparisonTMT

Group comparison for TMT data

Description

Group comparison for TMT data

Usage

```
MSstatsGroupComparisonTMT(fitted_models, contrast_matrix)
```

Arguments

fitted_models output of the MSstatsModerateTTest function

contrast_matrix

contrast matrix

Value

data.table

MSstatsModerateTTest *Moderate T statistic for group comparison*

Description

Moderate T statistic for group comparison

Usage

MSstatsModerateTTest(summarized, fitted_models, moderated)

Arguments

| | |
|---------------|--|
| summarized | protein-level data produced by the proteinSummarization function |
| fitted_models | output of the MSstatsFitComparisonModelsTMT function |
| moderated | if TRUE, moderation will be performed |

Value

list

MSstatsNormalizeTMT *Normalization for TMT data*

Description

Normalization for TMT data

Usage

MSstatsNormalizeTMT(input, type, normalize)

Arguments

| | |
|-----------|--|
| input | data.table |
| type | "peptides" for peptide normalization between channel and run, "proteins" for protein normalization |
| normalize | logical, if TRUE, data will be normalized |

Value

data.table

MSstatsPrepareForGroupComparisonTMT

Prepare output of proteinSummarization for group comparison

Description

Prepare output of proteinSummarization for group comparison

Usage

```
MSstatsPrepareForGroupComparisonTMT(  
  input,  
  remove_norm_channel,  
  remove_empty_channel  
)
```

Arguments

input output of proteinSummarization
remove_norm_channel if TRUE, "Norm" channel will be removed
remove_empty_channel if TRUE, empty channel will be removed

Value

data.table

MSstatsPrepareForSummarizationTMT

Prepare output of MSstatsTMT converters for protein-level summarization

Description

Prepare output of MSstatsTMT converters for protein-level summarization

Usage

```
MSstatsPrepareForSummarizationTMT(  
  data,  
  method,  
  global_norm,  
  reference_norm,  
  remove_norm_channel,
```

```

    remove_empty_channel,
    MBimpute,
    maxQuantileforCensored
  )

```

Arguments

| | |
|------------------------|---|
| data | Name of the output of PDtoMSstatsTMTFormat function or peptide-level quantified data from other tools. It should have columns ProteinName, PeptideSequence, Charge, PSM, Mixture, TechRepMixture, Run, Channel, Condition, BioReplicate, Intensity |
| method | Four different summarization methods to protein-level can be performed : "msstats"(default), "MedianPolish", "Median", "LogSum". |
| global_norm | Global median normalization on peptide level data (equalizing the medians across all the channels and MS runs). Default is TRUE. It will be performed before protein-level summarization. |
| reference_norm | Reference channel based normalization between MS runs on protein level data. TRUE(default) needs at least one reference channel in each MS run, annotated by 'Norm' in Condition column. It will be performed after protein-level summarization. FALSE will not perform this normalization step. If data only has one run, then reference_norm=FALSE. |
| remove_norm_channel | TRUE(default) removes 'Norm' channels from protein level data. |
| remove_empty_channel | TRUE(default) removes 'Empty' channels from protein level data. |
| MBimpute | only for method="msstats". TRUE (default) imputes missing values by Accelerated failure model. FALSE uses minimum value to impute the missing value for each peptide precursor ion. |
| maxQuantileforCensored | We assume missing values are censored. maxQuantileforCensored is Maximum quantile for deciding censored missing value, for instance, 0.999. Default is Null. |

Value

data.table

MSstatsSummarizationOutputTMT

Combine feature-level and protein-level data into single output

Description

Combine feature-level and protein-level data into single output

Usage

```
MSstatsSummarizationOutputTMT(
  summarized,
  processed,
  remove_empty_channel,
  remove_norm_channel
)
```

Arguments

summarized output of the getSummarizedTMT function
 processed output of the getProcessedTMT function
 remove_empty_channel
 TRUE(default) removes 'Empty' channels from protein level data.
 remove_norm_channel
 TRUE(default) removes 'Norm' channels from protein level data.

Value

list that consists of two data.frames with feature-level and protein-level data

MSstatsSummarizeTMT *Protein summarization for TMT data*

Description

Protein summarization for TMT data

Usage

```
MSstatsSummarizeTMT(
  input,
  method,
  impute,
  max_quantile_censored = NULL,
  log_file_path = NULL
)
```

Arguments

input data.table with TM quant data
 method Four different summarization methods to protein-level can be performed : "msstats"(default), "MedianPolish", "Median", "LogSum".
 impute only for method="msstats". TRUE (default) imputes missing values by Accelerated failure model. FALSE uses minimum value to impute the missing value for each peptide precursor ion.

max_quantile_censored
 We assume missing values are censored. maxQuantileforCensored is Maximum quantile for deciding censored missing value, for instance, 0.999. Default is Null.

log_file_path path to a MSstats log file

Value

data.table

MSstatsTestSingleProteinTMT
Hypothesis tests for a single protein in TMT data

Description

Hypothesis tests for a single protein in TMT data

Usage

```
MSstatsTestSingleProteinTMT(fitted_model, contrast_matrix)
```

Arguments

fitted_model single element of the MSstatsModerateTTest output

contrast_matrix
 contrast matrix

Value

list

MSstatsTMT *MSstatsTMT: A package for protein significance analysis in shotgun mass spectrometry-based proteomic experiments with tandem mass tag (TMT) labeling*

Description

A set of tools for detecting differentially abundant peptides and proteins in shotgun mass spectrometry-based proteomic experiments with tandem mass tag (TMT) labeling.

functions

- [PDtoMSstatsTMTFormat](#) : generates MSstatsTMT required input format for Proteome discoverer output.
- [MaxQtoMSstatsTMTFormat](#) : generates MSstatsTMT required input format for MaxQuant output.
- [SpectroMinetoMSstatsTMTFormat](#) : generates MSstatsTMT required input format for SpectroMine output.
- [OpenMStoMSstatsTMTFormat](#) : generates MSstatsTMT required input format for OpenMS output.
- [proteinSummarization](#) : summarizes PSM level quantification to protein level quantification.
- [dataProcessPlotsTMT](#) : visualizes for explanatory data analysis.
- [groupComparisonTMT](#) : tests for significant changes in protein abundance across conditions.

OpenMStoMSstatsTMTFormat

Generate MSstatsTMT required input format for OpenMS output

Description

Generate MSstatsTMT required input format for OpenMS output

Usage

```
OpenMStoMSstatsTMTFormat(
  input,
  useUniquePeptide = TRUE,
  rmPSM_withfewMea_withinRun = TRUE,
  rmProtein_with1Feature = FALSE,
  summaryforMultiplePSMs = sum,
  use_log_file = TRUE,
  append = FALSE,
  verbose = TRUE,
  log_file_path = NULL,
  ...
)
```

Arguments

| | |
|------------------|--|
| input | MSstatsTMT report from OpenMS |
| useUniquePeptide | TRUE(default) removes peptides that are assigned for more than one proteins. We assume to use unique peptide for each protein. |

| | |
|----------------------------|---|
| rmPSM_withfewMea_withinRun | TRUE (default) will remove the features that have 1 or 2 measurements within each Run. |
| rmProtein_with1Feature | TRUE will remove the proteins which have only 1 peptide and charge. Default is FALSE. |
| summaryforMultiplePSMs | sum(default) or max - when there are multiple measurements for certain feature in certain run, select the feature with the largest summation or maximal value. |
| use_log_file | logical. If TRUE, information about data processing will be saved to a file. |
| append | logical. If TRUE, information about data processing will be added to an existing log file. |
| verbose | logical. If TRUE, information about data processing will be printed to the console. |
| log_file_path | character. Path to a file to which information about data processing will be saved. If not provided, such a file will be created automatically. If 'append = TRUE', has to be a valid path to a file. |
| ... | additional parameters to 'data.table::fread'. |

Value

'data.frame' of class 'MSstatsTMT'.

Examples

```
head(raw.om)
input.om <- OpenMStoMSstatsTMTFormat(raw.om)
head(input.om)
```

PDtoMSstatsTMTFormat *Convert Proteome Discoverer output to MSstatsTMT format.*

Description

Convert Proteome Discoverer output to MSstatsTMT format.

Usage

```
PDtoMSstatsTMTFormat(
  input,
  annotation,
  which.proteinid = "Protein.Accessions",
  useNumProteinsColumn = TRUE,
  useUniquePeptide = TRUE,
  rmPSM_withfewMea_withinRun = TRUE,
```

```

    rmProtein_with1Feature = FALSE,
    summaryforMultipleRows = sum,
    use_log_file = TRUE,
    append = FALSE,
    verbose = TRUE,
    log_file_path = NULL,
    ...
)

```

Arguments

| | |
|----------------------------|---|
| input | PD report or a path to it. |
| annotation | annotation with Run, Fraction, TechRepMixture, Mixture, Channel, BioReplicate, Condition columns or a path to file. Refer to the example 'annotation' for the meaning of each column. |
| which.proteinid | Use 'Protein.Accessions' (default) column for protein name. 'Master.Protein.Accessions' can be used instead to get the protein name with single protein. |
| useNumProteinsColumn | logical, TRUE (default) remove shared peptides by information of # Proteins column in PSM sheet. |
| useUniquePeptide | logical, if TRUE (default) removes peptides that are assigned for more than one proteins. We assume to use unique peptide for each protein. |
| rmPSM_withfewMea_withinRun | TRUE (default) will remove the features that have 1 or 2 measurements within each Run. |
| rmProtein_with1Feature | TRUE will remove the proteins which have only 1 peptide and charge. Default is FALSE. |
| summaryforMultipleRows | sum (default) or max - when there are multiple measurements for certain feature in certain run, select the feature with the largest summation or maximal value. |
| use_log_file | logical. If TRUE, information about data processing will be saved to a file. |
| append | logical. If TRUE, information about data processing will be added to an existing log file. |
| verbose | logical. If TRUE, information about data processing will be printed to the console. |
| log_file_path | character. Path to a file to which information about data processing will be saved. If not provided, such a file will be created automatically. If 'append = TRUE', has to be a valid path to a file. |
| ... | additional parameters to 'data.table::fread'. |

Value

'data.frame' of class 'MSstatsTMT'

Examples

```
head(raw.pd)
head(annotation.pd)
input.pd <- PDtoMSstatsTMTFormat(raw.pd, annotation.pd)
head(input.pd)
```

PhilosophertoMSstatsTMTFormat

Convert Philosopher (Fragpipe) output to MSstatsTMT format.

Description

Convert Philosopher (Fragpipe) output to MSstatsTMT format.

Usage

```
PhilosophertoMSstatsTMTFormat(
  input,
  annotation,
  protein_id_col = "Protein",
  peptide_id_col = "Peptide.Sequence",
  Purity_cutoff = 0.6,
  PeptideProphet_prob_cutoff = 0.7,
  useUniquePeptide = TRUE,
  rmPSM_withfewMea_withinRun = TRUE,
  rmPeptide_OxidationM = TRUE,
  rmProtein_with1Feature = FALSE,
  summaryforMultipleRows = sum,
  use_log_file = TRUE,
  append = FALSE,
  verbose = TRUE,
  log_file_path = NULL,
  ...
)
```

Arguments

| | |
|----------------|--|
| input | data.frame of 'msstats.csv' file produced by Philosopher |
| annotation | annotation with Run, Fraction, TechRepMixture, Mixture, Channel, BioReplicate, Condition columns or a path to file. Refer to the example 'annotation' for the meaning of each column. Channel column should be consistent with the channel columns (Ignore the prefix "Channel ") in msstats.csv file. Run column should be consistent with the Spectrum.File columns in msstats.csv file. |
| protein_id_col | Use 'Protein'(default) column for protein name. 'Master.Protein.Accessions' can be used instead to get the protein ID with single protein. |

| | |
|----------------------------|---|
| peptide_id_col | Use 'Peptide.Sequence'(default) column for peptide sequence. 'Modified.Peptide.Sequence' can be used instead to get the modified peptide sequence. |
| Purity_cutoff | Cutoff for purity. Default is 0.6 |
| PeptideProphet_prob_cutoff | Cutoff for the peptide identification probability. Default is 0.7. The probability is confidence score determined by PeptideProphet and higher values indicate greater confidence. |
| useUniquePeptide | logical, if TRUE (default) removes peptides that are assigned for more than one proteins. We assume to use unique peptide for each proteins. |
| rmPSM_withfewMea_withinRun | TRUE (default) will remove the features that have 1 or 2 measurements within each Run. |
| rmPeptide_OxidationM | TRUE (default) will remove the peptides including oxidation (M) sequence. |
| rmProtein_with1Feature | TRUE will remove the proteins which have only 1 peptide and charge. Default is FALSE. |
| summaryforMultipleRows | sum (default) or max - when there are multiple measurements for certain feature in certain run, select the feature with the largest summation or maximal value. |
| use_log_file | logical. If TRUE, information about data processing will be saved to a file. |
| append | logical. If TRUE, information about data processing will be added to an existing log file. |
| verbose | logical. If TRUE, information about data processing will be printed to the console. |
| log_file_path | character. Path to a file to which information about data processing will be saved. If not provided, such a file will be created automatically. If 'append = TRUE', has to be a valid path to a file. |
| ... | additional parameters to 'data.table::fread'. |

Value

'data.frame' of class 'MSstatsTMT'

| | |
|---------------|--|
| proteinGroups | <i>Example of proteinGroups file from MaxQuant for TMT-10plex experiments.</i> |
|---------------|--|

Description

Example of proteinGroup.txt file from MaxQuant, which is identified protein group information file. It is the input for MaxQtoMSstatsTMTFormat function, with evidence.txt and annotation file. It includes identified protein groups for 10 proteins among 15 MS runs with TMT10. The important variables are as follows:

Usage

```
proteinGroups
```

Format

A data frame with 1075 rows and 105 variables.

Details

- id
- Protein.IDs
- Only.identified.by.site
- Potential.contaminant
- Reverse

Examples

```
head(proteinGroups)
```

proteinSummarization *Summarizing peptide level quantification to protein level quantification*

Description

We assume missing values are censored and then impute the missing values. Protein-level summarization from peptide level quantification are performed. After all, global median normalization on peptide level data and normalization between MS runs using reference channels will be implemented.

Usage

```
proteinSummarization(  
  data,  
  method = "msstats",  
  global_norm = TRUE,  
  reference_norm = TRUE,  
  remove_norm_channel = TRUE,  
  remove_empty_channel = TRUE,  
  MBimpute = TRUE,  
  maxQuantileforCensored = NULL,  
  use_log_file = TRUE,  
  append = FALSE,  
  verbose = TRUE,  
  log_file_path = NULL,  
  msstats_log_path = NULL  
)
```

Arguments

| | |
|------------------------|---|
| data | Name of the output of PDtoMSstatsTMTFormat function or peptide-level quantified data from other tools. It should have columns ProteinName, PeptideSequence, Charge, PSM, Mixture, TechRepMixture, Run, Channel, Condition, BioReplicate, Intensity |
| method | Four different summarization methods to protein-level can be performed : "msstats"(default), "MedianPolish", "Median", "LogSum". |
| global_norm | Global median normalization on peptide level data (equalizing the medians across all the channels and MS runs). Default is TRUE. It will be performed before protein-level summarization. |
| reference_norm | Reference channel based normalization between MS runs on protein level data. TRUE(default) needs at least one reference channel in each MS run, annotated by 'Norm' in Condition column. It will be performed after protein-level summarization. FALSE will not perform this normalization step. If data only has one run, then reference_norm=FALSE. |
| remove_norm_channel | TRUE(default) removes 'Norm' channels from protein level data. |
| remove_empty_channel | TRUE(default) removes 'Empty' channels from protein level data. |
| MBimpute | only for method="msstats". TRUE (default) imputes missing values by Accelerated failure model. FALSE uses minimum value to impute the missing value for each peptide precursor ion. |
| maxQuantileforCensored | We assume missing values are censored. maxQuantileforCensored is Maximum quantile for deciding censored missing value, for instance, 0.999. Default is Null. |
| use_log_file | logical. If TRUE, information about data processing will be saved to a file. |
| append | logical. If TRUE, information about data processing will be added to an existing log file. |
| verbose | logical. If TRUE, information about data processing will be printed to the console. |
| log_file_path | character. Path to a file to which information about data processing will be saved. If not provided, such a file will be created automatically. If 'append = TRUE', has to be a valid path to a file. |
| msstats_log_path | path to a MSstats log file |

Value

list that consists of two data.frames with feature-level (FeatureLevelData) and protein-level data (ProteinLevelData)

Examples

```
data(input.pd)
quant.pd.msstats <- proteinSummarization(input.pd,
                                         method = "msstats",
                                         global_norm = TRUE,
                                         reference_norm = TRUE)
head(quant.pd.msstats$ProteinLevelData)
```

quant.pd.msstats *Example of output from proteinSummarization function*

Description

It is made from [input.pd](#). It is the output of proteinSummarization function. It is a list that consists of two data.frames with feature-level (FeatureLevelData) and protein-level data (ProteinLevelData). ProteinLevelData should include the required columns as below.

Usage

```
quant.pd.msstats
```

Format

A data frame with 100 rows and 8 variables.

Details

- Run : MS run ID
- Protein : Protein ID
- Abundance: Protein-level summarized abundance
- Channel : Labeling information (126, ... 131)
- Condition : Condition (ex. Healthy, Cancer, Time0)
- BioReplicate : Unique ID for biological subject.
- TechRepMixture : Unique ID for technical replicate of one TMT mixture.
- Mixture : Unique ID for TMT mixture.

Examples

```
head(quant.pd.msstats$ProteinLevelData)
```

`raw.mine`*Example of output from SpectroMine for TMT-6plex experiments.*

Description

Example of SpectroMine PSM sheet. It is the output of SpectroMine and the input for SpectroMine-toMSstatsTMTFormat function, with annotation file. Annotation file should be made by users. It includes peak intensities for 10 proteins among 12 MS runs with TMT-6plex. The important variables are as follows:

Usage

`raw.mine`

Format

A data frame with 170 rows and 28 variables.

Details

- PG.ProteinAccessions
- P.MoleculeID
- PP.Charge
- R.FileName
- PG.QValue
- PSM.Qvalue
- Channels : PSM.TMT6_126..Raw., ..., PSM.TMT6_131..Raw.

Examples

`head(raw.mine)`

`raw.om`*Example of MSstatsTMT report from OpenMS for TMT-10plex experiments.*

Description

Example of MSstatsTMT PSM sheet from MaxQuant. It is the input for OpenMStoMSstatsTMT-Format function. It includes peak intensities for 10 proteins among 27 MS runs from three TMT10 mixtures. The important variables are as follows:

Usage

```
raw.om
```

Format

A data frame with 860 rows and 13 variables.

Details

- RetentionTime
- ProteinName
- PeptideSequence
- Charge
- Channel
- Condition
- BioReplicate
- Run
- Mixture
- TechRepMixture
- Fraction
- Intensity
- Reference

Examples

```
head(raw.om)
```

```
raw.pd
```

Example of output from Proteome Discoverer 2.2 for TMT-10plex experiments.

Description

Example of Proteome discover PSM sheet. It is the input for PDtoMSstatsTMTFormat function, with annotation file. Annotation file should be made by users. It includes peak intensities for 10 proteins among 15 MS runs with TMT-10plex. The variables are as follows:

Usage

```
raw.pd
```

Format

A data frame with 2858 rows and 50 variables.

Details

- Master.Protein.Accessions
- Protein.Accessions
- Annotated.Sequence
- Charge
- Ions.Score
- Spectrum.File
- Quan.Info
- Channels : 126, ..., 131

Examples

```
head(raw.pd)
```

SpectroMinetoMSstatsTMTFormat

Import data from SpectroMine

Description

Import data from SpectroMine

Usage

```
SpectroMinetoMSstatsTMTFormat(  
  input,  
  annotation,  
  filter_with_Qvalue = TRUE,  
  qvalue_cutoff = 0.01,  
  useUniquePeptide = TRUE,  
  rmPSM_withfewMea_withinRun = TRUE,  
  rmProtein_with1Feature = FALSE,  
  summaryforMultipleRows = sum,  
  use_log_file = TRUE,  
  append = FALSE,  
  verbose = TRUE,  
  log_file_path = NULL,  
  ...  
)
```

Arguments

| | |
|---|---|
| <code>input</code> | data name of SpectroMine PSM output. Read PSM sheet. |
| <code>annotation</code> | data frame which contains column Run, Fraction, TechRepMixture, Mixture, Channel, BioReplicate, Condition. Refer to the example 'annotation.mine' for the meaning of each column. |
| <code>filter_with_Qvalue</code> | TRUE(default) will filter out the intensities that have greater than <code>qvalue_cutoff</code> in EG.Qvalue column. Those intensities will be replaced with NA and will be considered as censored missing values for imputation purpose. |
| <code>qvalue_cutoff</code> | Cutoff for EG.Qvalue. default is 0.01. |
| <code>useUniquePeptide</code> | TRUE(default) removes peptides that are assigned for more than one proteins. We assume to use unique peptide for each protein. |
| <code>rmPSM_withfewMea_withinRun</code> | TRUE (default) will remove the features that have 1 or 2 measurements within each Run. |
| <code>rmProtein_with1Feature</code> | TRUE will remove the proteins which have only 1 peptide and charge. Default is FALSE. |
| <code>summaryforMultipleRows</code> | sum(default) or max - when there are multiple measurements for certain feature in certain run, select the feature with the largest summation or maximal value. |
| <code>use_log_file</code> | logical. If TRUE, information about data processing will be saved to a file. |
| <code>append</code> | logical. If TRUE, information about data processing will be added to an existing log file. |
| <code>verbose</code> | logical. If TRUE, information about data processing will be printed to the console. |
| <code>log_file_path</code> | character. Path to a file to which information about data processing will be saved. If not provided, such a file will be created automatically. If 'append = TRUE', has to be a valid path to a file. |
| <code>...</code> | additional parameters to 'data.table::fread'. |

Value

'data.frame' of class 'MSstatsTMT'

Examples

```
head(raw.mine)
head(annotation.mine)
input.mine <- SpectroMineMSstatsTMTFormat(raw.mine, annotation.mine)
head(input.mine)
```

`test.pairwise`*Example of output from groupComparisonTMT function*

Description

It is the output of groupComparisonTMT function, which is made from [quant.pd.msstats](#). It is a list that consists of the following elements: (1) ComparisonResult: statistical testing results; (2) FittedModel: the fitted linear models ComparisonResult should include the columns as below.

Usage

```
test.pairwise
```

Format

A data frame with 60 rows and 7 variables.

Details

- Protein : Protein ID
- Label: Label of the pairwise comparison or contrast
- log2FC: Log2 fold change
- SE: Standard error of the comparison of contrast results
- DF: Degree of freedom
- pvalue: Value of p statistic of the test
- adj.pvalue: adjusted p value
- issue: used for indicating the reason why a comparison is not testable. NA means the comparison is testable. 'oneConditionMissing' means the protein has no measurements in one condition of the comparison. Furthermore, when 'issue = oneConditionMissing', 'log2FC = Inf' means the negative condition (with coefficient -1 in the Label column) is missing and 'log2FC = -Inf' means the positive condition (with coefficient 1 in the Label column) is missing. 'completeMissing' means the protein has no measurements in all the conditions of the comparison. 'unfittableModel' means there is not enough measurements to fit the linear model. In other words, each condition has only one measurement.

Examples

```
head(test.pairwise$ComparisonResult)
```

Index

* datasets

- annotation.mine, 16
- annotation.mq, 17
- annotation.pd, 18
- evidence, 21
- input.pd, 24
- proteinGroups, 38
- raw.mine, 42
- raw.om, 42
- raw.pd, 43

* internal

- .checkContrastMatrix, 3
- .checkSummarizationParams, 3
- .countRunsWithNorm, 5
- .documentFunction, 5
- .getNormalizationAbundance, 6
- .getPhilosopherInput, 7
- .getRunsMedian, 7
- .handleSingleContrastTMT, 8
- .logSum, 8
- .logSummarizationParams, 9
- .makeContrastSingleTMT, 10
- .makeFactorColumnsTMT, 10
- .medianPolish, 11
- .normalizePeptides, 11
- .normalizeProteins, 12
- .prepareForSummarization, 12
- .removeRedundantChannels, 13
- .summarizeMSstats, 13
- .summarizeSimpleStat, 14
- .summarizeTMP, 15
- .summarizeTMT, 15
- getProcessedTMT, 22
- getSummarizedTMT, 22
- MSstatsComparisonModelSingleTMT, 27
- MSstatsFitComparisonModelsTMT, 27
- MSstatsGroupComparisonOutputTMT, 28

- MSstatsGroupComparisonTMT, 28
- MSstatsModerateTTest, 29
- MSstatsNormalizeTMT, 29
- MSstatsPrepareForGroupComparisonTMT, 30
- MSstatsPrepareForSummarizationTMT, 30
- MSstatsSummarizationOutputTMT, 31
- MSstatsSummarizeTMT, 32
- MSstatsTestSingleProteinTMT, 33
- quant.pd.msstats, 41
- test.pairwise, 46

- .checkContrastMatrix, 3
- .checkSummarizationParams, 3
- .countRunsWithNorm, 5
- .documentFunction, 5
- .getNormalizationAbundance, 6
- .getPhilosopherInput, 7
- .getRunsMedian, 7
- .handleSingleContrastTMT, 8
- .logSum, 8
- .logSummarizationParams, 9
- .makeContrastSingleTMT, 10
- .makeFactorColumnsTMT, 10
- .medianPolish, 11
- .normalizePeptides, 11
- .normalizeProteins, 12
- .prepareForSummarization, 12
- .removeRedundantChannels, 13
- .summarizeMSstats, 13
- .summarizeSimpleStat, 14
- .summarizeTMP, 15
- .summarizeTMT, 15

- annotation.mine, 16
- annotation.mq, 17
- annotation.pd, 18, 24

- dataProcessPlotsTMT, 19, 34

evidence, [21](#)

getProcessedTMT, [22](#)
getSummarizedTMT, [22](#)
groupComparisonTMT, [23](#), [34](#)

input.pd, [24](#), [41](#)

MaxQtoMSstatsTMTFormat, [25](#), [34](#)
MSstatsComparisonModelSingleTMT, [27](#)
MSstatsFitComparisonModelsTMT, [27](#)
MSstatsGroupComparisonOutputTMT, [28](#)
MSstatsGroupComparisonTMT, [28](#)
MSstatsModerateTTest, [29](#)
MSstatsNormalizeTMT, [29](#)
MSstatsPrepareForGroupComparisonTMT,
[30](#)
MSstatsPrepareForSummarizationTMT, [30](#)
MSstatsSummarizationOutputTMT, [31](#)
MSstatsSummarizeTMT, [32](#)
MSstatsTestSingleProteinTMT, [33](#)
MSstatsTMT, [33](#)

OpenMStoMSstatsTMTFormat, [34](#), [34](#)

PDtoMSstatsTMTFormat, [34](#), [35](#)
PhilosophertoMSstatsTMTFormat, [37](#)
proteinGroups, [38](#)
proteinSummarization, [19](#), [23](#), [34](#), [39](#)

quant.pd.msstats, [41](#), [46](#)

raw.mine, [42](#)
raw.om, [42](#)
raw.pd, [24](#), [43](#)

SpectroMinetoMSstatsTMTFormat, [34](#), [44](#)

test.pairwise, [46](#)