Package 'synapsis'

May 11, 2024

Type Package

Title An R package to automate the analysis of double-strand break repair during meiosis

Version 1.10.0

Description Synapsis is a Bioconductor software package for automated (unbiased and reproducible) analysis of meiotic immunofluorescence datasets. The primary functions of the software can i) identify cells in meiotic prophase that are labelled by a synaptonemal complex axis or central element protein, ii) isolate individual synaptonemal complexes and measure their physical length, iii) quantify foci and co-localise them with synaptonemal complexes, iv) measure interference between synaptonemal complex-associated foci. The software has applications that extend to multiple species and to the analysis of other proteins that label meiotic prophase chromosomes. The software converts meiotic immunofluorescence images into R data frames that are compatible with machine learning methods. Given a set of microscopy images of meiotic spread slides, synapsis crops images around individual single cells, counts colocalising foci on strands on a per cell basis, and measures the distance between foci on any given strand.

```
biocViews Software, SingleCell

Depends R (>= 4.1)

Imports EBImage, stats, utils, graphics

License MIT + file LICENSE

Encoding UTF-8

RoxygenNote 7.1.1

VignetteBuilder knitr

Suggests knitr, rmarkdown, testthat (>= 3.0.0), ggplot2, tidyverse, BiocStyle

Config/testthat/edition 3

git_url https://git.bioconductor.org/packages/synapsis

git_branch RELEASE_3_19

git_last_commit_2c703bf

git_last_commit_date 2024-04-30

Repository Bioconductor 3.19
```

Date/Publication 2024-05-10

```
Author Lucy McNeill [aut, cre, cph] (<a href="https://orcid.org/0000-0003-1752-4882">https://orcid.org/0000-0003-1752-4882</a>), Wayne Crismani [rev, ctb] (<a href="https://orcid.org/0000-0003-0143-8293">https://orcid.org/0000-0003-0143-8293</a>)
```

Maintainer Lucy McNeill < luc.mcneill@gmail.com>

Contents

```
20
Index
   23
annotate_foci_counting
 annotate_foci_counting
```

Description

Contains all plotting routines for count foci annotation

```
annotate_foci_counting(
  img_file,
  cell_count,
  img_orig,
  img_orig_foci,
  artificial_amp_factor,
  strands,
  coincident_foci,
  foci_label,
  alone_foci,
```

```
percent_px,
  foci_per_cell
)
```

Arguments

```
img_file
                 cell's file name
cell_count
                 unique cell counter
img_orig
                  original strand crop
img_orig_foci
                 cropped foci channel
artificial_amp_factor
                  amplification factor
                 black white mask of strand channel
strands
coincident_foci
                  mask of overlap between strand and foci channel
foci_label
                 black and white mask of foci channel
alone_foci
                  estimated number of foci that are NOT on a strand.
percent_px
                  percentage of foci mask that coincides with strand channel small number indi-
                 cates potentially problematic image.
foci_per_cell
                 number of foci counted per cell
```

Value

displays key steps from raw image to coincident foci count

```
annotate_foci_counting_adjusted

annotate_foci_counting_adjusted
```

Description

Contains all plotting routines for count foci annotation

```
annotate_foci_counting_adjusted(
  img_file,
  cell_count,
  img_orig,
  img_orig_foci,
  artificial_amp_factor,
  strands,
  coincident_foci,
  foci_label,
  alone_foci,
```

4 append_data_frame

```
percent_px,
  foci_per_cell
)
```

Arguments

img_file cell's file name cell_count unique cell counter img_orig original strand crop img_orig_foci cropped foci channel artificial_amp_factor amplification factor strands black white mask of strand channel coincident_foci mask of overlap between strand and foci channel foci_label black and white mask of foci channel alone_foci estimated number of foci that are NOT on a strand. percent_px percentage of foci mask that coincides with strand channel small number indicates potentially problematic image. number of foci counted per cell foci_per_cell

Value

displays key steps from raw image to coincident foci count

```
append_data_frame
append_data_frame
```

Description

applies new row to data frame

```
append_data_frame(
  WT_str,
  KO_str,
  WT_out,
  KO_out,
  img_file,
  foci_areas,
  df_cells,
  cell_count,
  stage,
  foci_per_cell,
```

auto_crop_fast 5

```
image_mat,
percent_px,
alone_foci,
discrepant_category,
C1
```

Arguments

WT_str string in filename corresponding to wildtype genotype. Defaults to ++.

KO_str string in filename corresponding to knockout genotype. Defaults to -.

WT_out string in output csv in genotype column, for knockout. Defaults to +/+.

KO_out string in output csv in genotype column, for knockout. Defaults to -/-.

img_file cell's file name

foci_areas pixel area of each foci
df_cells current data frame
cell_count unique cell counter

stage meiosis stage of interest. Currently count_foci determines this with threshold-

ing/ object properties in the synaptonemal complex channel by previously calling the get_pachytene function. Note that if using this option, the count_foci function requires that the input directory contains a folder called "pachytene" with

the crops in it.

foci_per_cell foci count for cell

image_mat matrix with all pixel values above zero

percent_px percentage of foci mask that coincides with strand channel small number indi-

cates potentially problematic image.

alone_foci estimated number of foci that are NOT on a strand.

discrepant_category

estimated number of foci that are NOT on a strand.

C1 criteria

Value

data frame with new row

auto_crop_fast

Description

crop an image around each viable cell candidate.

6 auto_crop_fast

Usage

```
auto_crop_fast(
  img_path,
 max_cell_area = 20000,
 min_cell_area = 7000,
 mean_pix = 0.08,
 annotation = "off",
 blob_factor = 15,
 bg_blob_factor = 10,
  offset = 0.2,
  final_blob_amp = 10,
  test_amount = 0,
  brush_size_blob = 51,
  sigma_blob = 15,
  channel3_string = "DAPI",
  channel2_string = "SYCP3"
  channel1_string = "MLH3",
  file_ext = "jpeg",
  third_channel = "off",
  cell_aspect_ratio = 2,
  strand_amp = 2,
  path_out = img_path,
  resize_1 = 720,
  crowded_cells = "FALSE",
 watershed_radius = 50,
 watershed_tol = 0.2,
  cropping_factor = 1.3
)
```

Arguments

img_path path containing image data to analyse max_cell_area Maximum pixel area of a cell candidate min_cell_area Minimum pixel area of a cell candidate Mean pixel intensity of cell crop (in SYCP3 channel) for normalisation mean_pix annotation Choice to output pipeline choices (recommended to knit) blob_factor Contrast factor to multiply original image by before smoothing/smudging bg_blob_factor Contrast factor to multiply original image by to take background. Used prior to thresholding. Pixel value offset from bg_blob_factor. Used in thresholding to make blob mask. offset Contrast factor to multiply smoothed/smudged image. Used in thresholding to final_blob_amp make blob mask. Optional number of first N images you want to run function on. For troutest_amount bleshooting/testing/variable calibration purposes. brush_size_blob

Brush size for smudging the synaptonemal complex channel to make blobs

auto_crop_fast 7

sigma_blob Sigma in Gaussian brush for smudging the synaptonemal complex channel to make blobs

channel3_string

Optional. String appended to the files showing the channel illuminating cell structures. Defaults to DAPI, if third channel == "on".

channel2_string

String appended to the files showing the channel illuminating synaptonemal complexes. Defaults to SYCP3

channel1_string

String appended to the files showing the channel illuminating foci. Defaults to

MLH3

file_ext file extension of your images e.g. tif jpeg or png.

third_channel Optional, defaults to "off". Set to "on" if you would also like crops of the third

channel.

cell_aspect_ratio

Maximum aspect ratio of blob to be defined as a cell

strand_amp multiplication of strand channel for get_blobs function.

path_out user specified output path. Defaults to img_path

resize_l length for resized image

crowded_cells TRUE or FALSE, defaults to FALSE. Set to TRUE if you have many cells in a

frame that almost touch

watershed_radius

Radius (ext variable) in watershed method used in strand channel. Defaults to 1

(small)

watershed_tol Intensity tolerance for watershed method. Defaults to 0.05.

cropping_factor

size of cropping window square, as factor of characteristic blob radius. Defaults to 1. May need to increase if using watershed.

Details

This function takes all images in a directory, and crops around individual cells according to the antibody that stains synaptonemal complexes e.g. SYCP3. First, it increases the brightness and smudges the image with a Gaussian brush, and creates a mask using thresholding (get_blobs). Then it deletes cell candidates in the mask deemed too large, too small, or too long (keep_cells). Using the computeFeatures functions from EBImage to locate centre and radius, the cropping area is determined and the original image cropped. These images are saved in either a user specified directory, or a crops folder at the location of the image files.

Value

cropped synaptonemal complex and foci channels around single cells, regardless of stage

Author(s)

Lucy McNeill

8 count_foci

Examples

```
demo_path = paste0(system.file("extdata",package = "synapsis"))
auto_crop_fast(demo_path, annotation = "on", max_cell_area = 30000,
min_cell_area = 7000, file_ext = "tif",crowded_cells = TRUE)
```

count_foci

count_foci

Description

Calculates coincident foci in synaptonemal complex and foci channel, per cell

```
count_foci(
  img_path,
  stage = "none",
  offset_px = 0.2,
  offset_factor = 2,
  brush\_size = 3,
  brush_sigma = 3,
  foci_norm = 0.01,
  annotation = "off"
  channel2_string = "SYCP3",
  channel1_string = "MLH3",
  file_ext = "jpeg",
  KO_str = "--",
  WT_str = "++"
 KO_out = "-/-"
 WT_out = "+/+",
 watershed_stop = "off",
 watershed_radius = 1,
 watershed_tol = 0.05,
  crowded_foci = TRUE,
  artificial_amp_factor = 1,
  strand_amp = 2,
 min_foci = -1,
  disc_size = 51,
 modify_problematic = "off",
  disc_size_foci = 5,
  C1 = 0.02,
  C2 = 0.46,
  C_weigh_foci_number = TRUE
)
```

count_foci 9

Arguments

img_path path containing crops to analyse

stage meiosis stage of interest. Currently count_foci determines this with threshold-

ing/ object properties in the synaptonemal complex channel by previously calling the get_pachytene function. Note that if using this option, the count_foci function requires that the input directory contains a folder called "pachytene" with

the crops in it.

offset_px Pixel value offset used in thresholding of synaptonemal complex channel

offset_factor Pixel value offset used in thresholding of foci channel

brush_size size of brush to smooth the foci channel. Should be small to avoid erasing foci.

brush_sigma sigma for Gaussian smooth of foci channel. Should be small to avoid erasing

foci.

foci_norm Mean intensity to normalise all foci channels to.

annotation Choice to output pipeline choices (recommended to knit)

channel2_string

String appended to the files showing the channel illuminating synaptonemal

complexes. Defaults to SYCP3

channel1_string

String appended to the files showing the channel illuminating foci. Defaults to

MLH3

file_ext file extension of your images e.g. tiff jpeg or png.

KO_str string in filename corresponding to knockout genotype. Defaults to -.

WT_str string in filename corresponding to wildtype genotype. Defaults to ++.

KO_out string in output csv in genotype column, for knockout. Defaults to -/-.

WT_out string in output csv in genotype column, for knockout. Defaults to +/+.

watershed_stop Stop default watershed method with "on"

watershed_radius

Radius (ext variable) in watershed method used in foci channel. Defaults to 1

(small)

watershed_tol Intensity tolerance for watershed method. Defaults to 0.05.

crowded_foci TRUE or FALSE, defaults to FALSE. Set to TRUE if you have foci > 100 or so.

artificial_amp_factor

Amplification of foci channel, for annotation only.

strand_amp multiplication of strand channel to make masks

min_foci minimum pixel area for a foci. Depends on your dpi etc. Defaults to 4

disc_size size of disc for local background calculation in synaptonemal complex channel

modify_problematic

option for synapsis to try and "save" images which have likely been counted incorrectly due to a number of reasons. Default settings are optimized for mouse

pachytene. Defaults to "off"

disc_size_foci size of disc for local background calculation in foci channel

Details

In this function, masks for the synaptonemal complex (SC) and foci channel are created from the saved crops of single/individual cells. These masks are computed using (optional) input parameters related to meiosis stage/ how well spread chromosomes are (for the former) and related to smoothing, thresholding and how "crowded" foci are for the latter. Finally, these two masks are multiplied, and the number of objects found with EBImage's computeFeatures are the colocalizing foci.

The file, cell number, foci count etc. are output as a data frame.

Value

data frame with foci count per cell

Author(s)

Lucy McNeill

Examples

Description

Creates mask for every individual cell candidate in mask

```
crop_single_object_fast(
  retained,
  OOI_final,
  counter_final,
  img_orig,
  img_orig_foci,
  img_orig_DAPI = "blank",
  file_sc,
```

```
file_foci,
  file_DAPI = "blank",
  cell_count,
 mean_pix,
  annotation,
  file_base,
  img_path,
  r_max,
  CX,
  су,
  channel3_string,
  channel2_string,
  channel1_string,
  file_ext,
  third_channel,
  path_out,
  img_orig_highres,
  resize_1,
  crowded_cells,
  cropping_factor
)
```

Arguments

retained Mask of cell candidates which meet size criteria. After smoothing/smudging

and thresholding.

OOI_final Objects of interest count. Total number of cell candidates in retained.

counter_final Counter for single cell we are focussing on. Remove all other cells where

counter_single not equal to counter_final.

img_orig description
img_orig_foci description
img_orig_DAPI description

file_sc filename of synaptonemal complex channel image

file_foci filename of foci channel image file_DAPI filename of DAPI channel image

cell_count counter for successful crops around cells

mean_pix Mean pixel intensity of cell crop (in SYCP3 channel) for normalisation

annotation Choice to output pipeline choices (recommended to knit)

file_base filename base common to all three channels i.e. without -MLH3.jpeg etc.

img_pathpath containing image data to analyser_maxmaximum radius of blob for cropping

cx centre of blob x cy centre of blob y

12 get_blobs

channel3_string

Optional. String appended to the files showing the channel illuminating cell structures. Defaults to DAPI, if third channel == "on".

channel2_string

String appended to the files showing the channel illuminating synaptonemal complexes. Defaults to SYCP3

channel1_string

String appended to the files showing the channel illuminating foci. Defaults to

MLH3

file_ext file extension of your images e.g. tif jpeg or png.

third_channel Optional, defaults to "off". Set to "on" if you would also like crops of the third

channel.

path_out user specified output path. Defaults to img_path

img_orig_highres

the original strand image with original resolution

resize_l length of square to resize original image to.

crowded_cells TRUE or FALSE, defaults to FALSE. Set to TRUE if you have many cells in a

frame that almost touch

cropping_factor

size of cropping window square, as factor of characteristic blob radius. Defaults to 1. May need to increase if using watershed.

Value

Crops around all candidates in both channels

get_blobs

get_blobs

Description

Makes mask of all objects bright enough to be classified as a cell cadidate

```
get_blobs(
  img_orig,
  blob_factor,
  bg_blob_factor,
  offset,
  final_blob_amp,
  brush_size_blob,
  sigma_blob,
  watershed_tol,
  watershed_radius,
  crowded_cells,
  annotation
)
```

get_C1 13

Arguments

img_orig Original image

blob_factor Contrast factor to multiply original image by before smoothing/smudging

bg_blob_factor Contrast factor to multiply original image by to take background. Used prior to

thresholding.

offset Pixel value offset from bg_blob_factor. Used in thresholding to make blob mask.

final_blob_amp Contrast factor to multiply smoothed/smudged image. Used in thresholding to

make blob mask.

brush_size_blob

Brush size for smudging the synaptonemal complex channel to make blobs

sigma_blob Sigma in Gaussian brush for smudging the synaptonemal complex channel to

make blobs

watershed_tol Intensity tolerance for watershed method. Defaults to 0.05.

watershed_radius

Radius (ext variable) in watershed method used in strand channel. Defaults to 1

(small)

crowded_cells TRUE or FALSE, defaults to FALSE. Set to TRUE if you have many cells in a

frame that almost touch

annotation Choice to output pipeline choices (recommended to knit) have many cells in a

frame that almost touch

Value

Mask with cell candidates

get_C1 get_C1

Description

calculates the statistic to compare to crisp_criteria, which determines whether the foci count will be reliable

Usage

```
get_C1(foci_areas, foci_per_cell, C_weigh_foci_number)
```

Arguments

foci_areas pixel area of each foci foci_per_cell foci count for cell C_weigh_foci_number

choose crispness criteria- defaults to TRUE to use C1 (weighing with number). Otherwise set to FALSE to use C2

14 get_coincident_foci

Value

statistic to comapre to crisp_criteria

Description

calculates the statistic to compare to crisp_criteria, which determines whether the foci count will be reliable

```
get_coincident_foci(
 offset_px,
  offset_factor,
 brush_size,
  brush_sigma,
  annotation,
 watershed_stop,
 watershed_radius,
 watershed_tol,
  crowded_foci,
  artificial_amp_factor,
  strand_amp,
  disc_size,
  disc_size_foci,
  img_file,
  cell_count,
  img_orig,
  img_orig_foci,
  stage,
 WT_str,
 KO_str,
 WT_out,
 KO_out,
 C1_search,
 discrepant_category,
 C1,
 C2,
 df_cells,
  C_weigh_foci_number
)
```

get_coincident_foci 15

Arguments

offset_px Pixel value offset used in thresholding of synaptonemal complex channel

offset_factor Pixel value offset used in thresholding of foci channel

brush_size size of brush to smooth the foci channel. Should be small to avoid erasing foci.

brush_sigma sigma for Gaussian smooth of foci channel. Should be small to avoid erasing

foci.

annotation Choice to output pipeline choices (recommended to knit)

watershed_stop Stop default watershed method with "on"

watershed_radius

Radius (ext variable) in watershed method used in foci channel. Defaults to 1

(small)

watershed_tol Intensity tolerance for watershed method. Defaults to 0.05.

crowded_foci TRUE or FALSE, defaults to FALSE. Set to TRUE if you have foci > 100 or so.

artificial_amp_factor

Amplification of foci channel, for annotation only.

strand_amp multiplication of strand channel to make masks

disc_size size of disc for local background calculation in synaptonemal complex channel

disc_size_foci size of disc for local background calculation in foci channel

img_file cell's file name
cell_count unique cell counter
img_orig original strand crop
img_orig_foci cropped foci channel

stage meiosis stage of interest. Currently count_foci determines this with threshold-

ing/object properties in the synaptonemal complex channel by previously calling the get_pachytene function. Note that if using this option, the count_foci function requires that the input directory contains a folder called "pachytene" with

the crops in it.

WT_str string in filename corresponding to wildtype genotype. Defaults to ++.

KO_str string in filename corresponding to knockout genotype. Defaults to -.

WT_out string in output csv in genotype column, for knockout. Defaults to +/+.

KO_out string in output csv in genotype column, for knockout. Defaults to -/-.

C1_search TRUE or FALSE whether the image is still being modified until it meets the

crispness criteria

discrepant_category

estimated number of foci that are NOT on a strand.

C1 Default crispness criteria = sd(foci_area)/(mean(foci_area)+1)

C2 Alternative crisp criteria.

df_cells current data frame

C_weigh_foci_number

choose crispness criteria- defaults to TRUE to use C1 (weighing with number).

Otherwise set to FALSE to use C2

get_foci_per_cell

Value

data frame with new row with most recent foci per cell appended

Description

creates mask for coincident foci

Usage

```
get_foci_per_cell(
   img_file,
   offset_px,
   stage,
   strands,
   watershed_stop,
   foci_label,
   annotation,
   cell_count,
   img_orig,
   img_orig_foci,
   artificial_amp_factor,
   coincident_foci
)
```

Arguments

img_file

_	
offset_px	Pixel value offset used in thresholding of synaptonemal complex channel

stage meiosis stage of interest. Currently count_foci determines this with thresholding/object properties in the synantonemal complex channel by previously calling

ing/ object properties in the synaptonemal complex channel by previously calling the get_pachytene function. Note that if using this option, the count_foci function requires that the input directory contains a folder called "pachytene" with

the crops in it.

cell's file name

strands black white mask of strand channel
watershed_stop Stop default watershed method with "on"
foci_label black and white mask of foci channel

annotation Choice to output pipeline choices (recommended to knit)

cell_count unique cell counter
img_orig original strand crop
img_orig_foci cropped foci channel

get_overlap_mask 17

Value

number of foci per cell

get_overlap_mask
get_overlap_mask

Description

creates mask for coincident foci

Usage

```
get_overlap_mask(
    strands,
    foci_label,
    watershed_stop,
    img_orig_foci,
    watershed_radius,
    watershed_tol
)
```

Arguments

strands black white mask of strand channel

foci_label black and white mask of foci channel

watershed_stop Stop default watershed method with "on"

img_orig_foci cropped foci channel

watershed_radius

Radius (ext variable) in watershed method used in foci channel. Defaults to 1

(small)

watershed_tol Intensity tolerance for watershed method. Defaults to 0.05.

Value

mask with coincident foci on strands

18 get_pachytene

get_pachytene

get_pachytene

Description

Identifies crops in pachytene

Usage

```
get_pachytene(
  img_path,
  species_num = 20,
 offset = 0.2,
  ecc_{thresh} = 0.85,
  area_thresh = 0.06,
  annotation = "off",
  channel2_string = "SYCP3",
  channel1_string = "MLH3",
  file_ext = "jpeg",
 KO_str = "--",
 WT_str = "++"
 KO_{out} = "-/-"
 WT_out = "+/+",
  path_out = img_path,
  artificial_amp_factor = 3,
  strand_amp = 2,
  resize_1 = 120
)
```

MLH3

Arguments

img_path path containing crops analyse species_num number of chromosomes in the species offset Pixel value offset used in therholding for the synaptonemal complex (SYCP3) channel The minimum average eccentricity of all objects in mask determined by comecc_thresh putefeatures, for a cell to be pachytene. The minimum ratio of pixels included in mask to total, for a cell to be classified area_thresh as pachytene. Choice to output pipeline choices (recommended to knit) annotation channel2_string String appended to the files showing the channel illuminating synaptonemal complexes. Defaults to SYCP3 channel1_string String appended to the files showing the channel illuminating foci. Defaults to keep_cells 19

file_ext	file extension of your images e.g. tiff jpeg or png.		
KO_str	string in filename corresponding to knockout genotype. Defaults to		
WT_str	string in filename corresponding to wildtype genotype. Defaults to ++		
KO_out	string in output csv in genotype column, for knockout. Defaults to -/		
WT_out	string in output csv in genotype column, for knockout. Defaults to +/+		
path_out	user specified output path. Defaults to img_path		
artificial_amp_factor Amplification of foci channel, for RGB output files. Deaults to 3.			
strand_amp	multiplication of strand channel.		
resize_l	length of resized square cell image.		

Details

This function takes the crops make by auto_crop fast, and determines the number of synaptonemal complex candidates by considering the local background and using EBImage functions. In general, very bright objects which contrast highly with the background will be classified as the same object. Dim objects will likely be classified as many different objects. If the number of objects is too high compared to the species number (species_num) then the cell is determined to not be in pachytene. Note that this function has been optimized for mouse cells which can be very well spread / separated.

Value

Pairs of foci and synaptonemal channel crops for pachytene

Author(s)

Lucy McNeill

Examples

```
demo_path = paste0(system.file("extdata",package = "synapsis"))
SYCP3_stats <- get_pachytene(demo_path,ecc_thresh = 0.8, area_thresh = 0.04, annotation = "on")</pre>
```

keep_cells		
	keep_cells	keep_cells

Description

Deletes objects in mask which are too small, large, oblong i.e. unlikely to be a cell

20 make_foci_mask

Usage

```
keep_cells(
  candidate,
  max_cell_area,
  min_cell_area,
  cell_aspect_ratio,
  crowded_cells,
  annotation
)
```

Arguments

Value

Mask of cell candidates which meet size criteria

Description

creates foci mask for foci channel crop

```
make_foci_mask(
   offset_factor,
   bg,
   crowded_foci,
   img_orig_foci,
   brush_size,
   brush_sigma,
   disc_size_foci
)
```

21 make_strand_mask

Arguments

offset_factor Pixel value offset used in thresholding of foci channel background value- currently just mean pixel value of whole image crowded_foci TRUE or FALSE, defaults to FALSE. Set to TRUE if you have foci > 100 or so. img_orig_foci cropped foci channel brush_size size of brush to smooth the foci channel. Should be small to avoid erasing foci. brush_sigma

sigma for Gaussian smooth of foci channel. Should be small to avoid erasing

disc_size_foci size of disc for local background calculation in foci channel

Value

foci mask

make_strand_mask

make_strand_mask

Description

creates strand mask for strand channel crop

Usage

```
make_strand_mask(
 offset_px,
  stage,
  img_orig,
  disc_size,
  brush_size,
  brush_sigma
)
```

Arguments

brush_sigma

offset_px Pixel value offset used in thresholding of synaptonemal complex channel meiosis stage of interest. Currently count_foci determines this with thresholdstage ing/object properties in the synaptonemal complex channel by previouly calling the get_pachytene function. Note that if using this option, the count_foci function requires that the input directory contains a folder called "pachytene" with the crops in it. img_orig original strand crop disc_size size of disc for local background calculation in synaptonemal complex channel brush_size size of brush to smooth the foci channel. Should be small to avoid erasing foci.

sigma for Gaussian smooth of foci channel. Should be small to avoid erasing

foci.

22 remove_XY

Value

strand mask

remove_XY

remove_XY

Description

applies new row to data frame

Usage

```
remove_XY(foci_label, foci_candidates, foci_areas)
```

Arguments

foci_label black and white mask of foci channel

 ${\tt foci_candidates}$

computeFeatures data frame of foci channel

foci_areas the areas of the foci objects

Value

mask with XY blob removed

Index

```
annotate\_foci\_counting, 2
{\tt annotate\_foci\_counting\_adjusted, 3}
append_data_frame, 4
auto_crop_fast, 5
count_foci, 8
{\tt crop\_single\_object\_fast, 10}
{\tt get\_blobs}, {\tt 12}
get_C1, 13
get_coincident_foci, 14
get_foci_per_cell, 16
get_overlap_mask, 17
get_pachytene, 18
keep_cells, 19
make\_foci\_mask, 20
{\tt make\_strand\_mask}, \textcolor{red}{21}
remove_XY, 22
```