

Enzyme	Rebase Motif	Example sequence
EcoRI	G^AATTC	<div>5'—ACGAACGAATTCGAGCGAGAAGG—3'</div> <div>3'—TGCTTGCTTAAGCTCGCTCTTCC—5'</div>
SmaI	CCC^GGG	<div>5'—ACGAACCCCGGGGAGCGAGAAGG—3'</div> <div>3'—TGCTTGGGGCCCTCGCTCTTCC—5'</div>
HgaI	GACGC(5/10)	<div>5'—ACGAACGACGCGGAGCGAGAAGG—3'</div> <div>3'—TGCTTGCTGCGCCTCGCTCTTCC—5'</div>
PfaI	G^GYRCC	<div>5'—ACGAACGGYRCCGAGCGAGAAGG—3'</div> <div>3'—TGCTTGCCRYGGCTCGCTCTTCC—5'</div>