Package 'DNAcopy'

May 17, 2024

Title DNA Copy Number Data Analysis **Version** 1.78.0 Author Venkatraman E. Seshan, Adam Olshen **Description** Implements the circular binary segmentation (CBS) algorithm to segment DNA copy number data and identify genomic regions with abnormal copy number. Maintainer Venkatraman E. Seshan <seshanv@mskcc.org> LazyData yes License GPL (>= 2) biocViews Microarray, CopyNumberVariation NeedsCompilation yes git_url https://git.bioconductor.org/packages/DNAcopy git_branch RELEASE_3_19 git_last_commit 51f8d46 git_last_commit_date 2024-04-30 **Repository** Bioconductor 3.19

Date/Publication 2024-05-17

Contents

CNA	2
coriell	3
cytoBand	4
DNAcopy	4
DNAcopy-internal	5
exon.segment	6
getbdry	7
glFrequency	7
plot.DNAcopy	8
plotSample	10
segment	12

segments.p		•	 •		•							•				•	•	•	15
segments.summary	•																		16
smooth.CNA																			17
subset.CNA																			18
subset.DNAcopy .																			19
zoomIntoRegion .																			20
-																			
																			22

Index

CNA

Create 'Copy Number Array' data object

Description

Creates a 'copy number array' data object used for DNA copy number analyses by programs such as circular binary segmentation (CBS).

Usage

genomdat	a vector or matrix of data from array-CGH, ROMA, or other copy number ex- periments. If it is a matrix the rows correspond to the markers and the columns to the samples.
chrom	the chromosomes (or other group identifier) from which the markers came. Vec- tor of length same as the number of rows of genomdat. If one wants the chro- mosomes to be ordered in the natural order, this variable should be numeric or ordered category.
maploc	the locations of marker on the genome. Vector of length same as the number of rows of genomdat. This has to be numeric.
data.type	logratio (aCGH, ROMA, etc.) or binary (LOH).
sampleid	sample identifier. If missing the samples are named by prefixing "Sample" to consecutive integers.
presorted	logical indicator telling if the data have already been sorted by chrom and maploc. Default is FALSE.
x	object returned by CNA
	arguments to be passed onto print command called within.

coriell

Details

Data that are NA, Inf, NaN will be removed on a per sample basis for "genomdat" and all samples for "chrom" and "maploc".

If the chrom variable has non-numeric values make it into an ordered variable to get them ordered correctly. E.g. for human genome use: chrom <- ordered(chrom, levels=c(1:22, "X", "Y")) to prepare the variable if chromosomes X and Y are present in your data.

Value

An object of class CNA. There is a print method that gives the number of samples and probes and the type of data.

Examples

coriell

Array CGH data set of Coriell cell lines

Description

These are two data array CGH studies sets of Corriel cell lines taken from the reference below.

Usage

data(coriell)

Format

A data frame containing five variables: first is clone name, second is clone chromosome, third is clone position, fourth (sample GM05296 from Web Table G) and fifth (sample GM13330 from Web Table H) are log2ratio for two cell lines.

Source

Supplementary Information from https://www.nature.com/articles/ng754

References

Snijders et al., Assembly of microarrays for genome-wide measurement of DNA copy number, *Nature Genetics*, 2001

cytoBand

Description

Cytogenic band data from the goldenPath repository

Usage

data(cytoBand)

Format

A data frame containing five variables: chromosome, start and end positions, band name and giesma stain.

Source

http://hgdownload.cse.ucsc.edu/goldenPath/hg19/database/cytoBand.txt.gz

DNAcopy

Results of segmenting a CNA data object

Description

The results of segmenting data from copy number array experiments from programs such as circular binary segmentation (CBS).

Usage

S3 method for class 'DNAcopy'
print(x, showSegRows=FALSE, ...)

Arguments

х	an object of class DNAcopy – output of segment.
showSegRows	option to show row numbers for the segment start and end. default is FALSE.
	arguments to be passed onto print command called within.

Details

An object of class DNAcopy. There is a print method that prints the results in a tabular format. Each row gives the sample, the chromosome, the start and end map locations, the number of markers and the mean of each segment.

DNAcopy-internal

Value

data	The original CNA object which was the input for segment
ID	sample identifier.
chrom	the chromosome within the sample.
loc.start	the starting map location of the segment
loc.end	the ending map location of the segment
num.mark	the number of markers in the segment
data.type	the segment mean.
call	the call that produced the object.

DNAcopy-internal

Internal DNAcopy functions

Description

Internal functions of package DNAcopy.

Usage

Details

These are not to be called directly by the user

exon.segment

Description

Compute the binary segmentation statistic, location and approximate p-value.

Usage

```
exon.segment(gene, eloc, edat, ngrid=100, tol=1e-6)
```

Arguments

gene names in the exon data
exon locations within gene
exon expressions within gene
number grid points for the integral
tolerance level for calculating nu

Details

The p-values are obtained by applying Siegmund's approximation for the maximal statistic from binary segmenting consecutive segments within a chromosome. These are one-sided test for an increase in expression.

Value

a matrix with three columns. The maximal statistic from binary segmentation, its location and the p-values for each gene.

Author(s)

Venkatraman E. Seshan

Examples

```
# test code on an easy data set
set.seed(25)
gene <- rep(c("A", "B"), c(30,20))
eloc <- c(1:30, 1:20)
edat <- matrix(rnorm(500), 50, 10)
# changes for gene1 in samples 3 & 7
edat[1:30, 3] <- edat[1:30, 3] + rep(0.9*0:1, c(17, 13))
edat[1:30, 7] <- edat[1:30, 7] + rep(1.1*0:1, c(21, 9))
# changes for gene2 in samples 4 & 7
edat[31:50, 4] <- edat[31:50, 4] + rep(1.1*0:1, c(8, 12))
edat[31:50, 7] <- edat[31:50, 7] + rep(1.2*0:1, c(13, 7))
exon.segment(gene, eloc, edat)
```

getbdry

Description

Function to compute the sequential boundary for early stopping.

Usage

getbdry(eta, nperm, max.ones, tol= 1e-2)

Arguments

eta	Type I error rate of the boundary.
nperm	Number of permutations for the reference distribution.
max.ones	maximum number of ones given by "floor(nperm*alpha)+1".
tol	tolerance level for the iterations.

Value

A vector integer values of length max.ones*(max.ones+1)/2 corresponding to the boundary for the number of ones from 1 to max.ones. The default boundary for nperm=10000, eta=0.05, alpha=0.01 is stored in the data object "default.DNAcopy.bdry". Use this function to get the boundary for your favorite values for the parameters "nperm, eta, alpha" and use it for the argument "sbdry" in the function "segment."

glFrequency

Additional summary measured for the segments

Description

This program computes the frequency of gains and losses for each probe as a function of level of mad.

Usage

glFrequency(xout, threshold=1)

xout	an object of class DNAcopy
threshold	threshold value to call gain or loss

Value

A segment is called a gain or loss if the segment mean is at least the threshold* mad distance away from the median copy number level. The output is a data frame with five columns which give the chromosome (chrom), genomic position (maploc), the number of samples with available data (pfreq), and the gain (gain) and loss (loss).

Author(s)

Venkatraman E. Seshan

plot.DNAcopy

Plot the data and results from segment of a CNA object

Description

Plots the data from a copy number array experiment (aCGH, ROMA etc.) along with the results of segmenting it into regions of equal copy numbers.

Usage

x	an object of class DNAcopy resulting from analyzing data from copy number array experiments.
plot.type	the type of plot.
xmaploc	logical flag to indicate that the X axis is the maploc position rather than the index. Since the segments are rearranged the plateau plot does not use maploc position.
altcol	logical flag to indicate if chromosomes should be plotted in alternating colors in the whole genome plot.
sbyc.layout	layout settings for the multifigure grid layout for the 'samplebychrom' type. It should be specified as a vector of two integers which are the number of rows and columns. The default values are chosen based on the number of chromosomes to produce a near square graph. For normal genome it is $4x6$ (24 chromosomes) plotted by rows.
cbys.layout	layout settings for the multifigure grid layout for the 'chrombysample' type. As above it should be specified as number of rows and columns and the default chosen based on the number of samples.

plot.DNAcopy

cbys.nchrom	the number of chromosomes per page in the layout. The default is 1.
include.means	logical flag to indicate whether segment means are to be drawn.
zeroline	logical flag to indicate whether a horizontal line at y=0 is to be drawn.
pt.pch	the plotting character used for plotting the log-ratio values (default is ".").
pt.cex	the size of plotting character used for the log-ratio values (default is 3).
pt.cols	the color list for the points. The colors alternate between chromosomes. If missing the point colors are black and green.
segcol	the color of the lines indicating the segment means. If missing the line color is set to be red.
zlcol	the color of the zeroline. If missing it is set to be grey.
ylim	this argument is present to override the default limits which is the range of sym- metrized log-ratios.
lwd	line weight of lines for segment mean and zeroline. If missing it is set to 3.
	other arguments which will be passed to plot commands.

Details

There are four possible plot types. For the type 'whole' the data are plotted for the entire genome. For the 'samplebychrom' type a graph with each chromosome (of a given sample) is drawn in a separate figure on a multi-figure grid. For the 'plateau' type the graph is drawn with the chromosome segments re-ordered by the segment means. For the 'chrombysample' type the samples for a given chromosome are drawn in a 4x6 multi-figure grid in multiples of 24. By default the segments means are drawn. For multisample data each sample or chromosome is drawn on a separate sheet. When invoked interactively the user is prompted before advancing to the next sample.

Examples

```
plot(segment.smoothed.CNA.object, plot.type="w")
#Plot each study by chromosome
plot(segment.smoothed.CNA.object, plot.type="s")
#Plot each chromosome across studies (6 per page)
plot(segment.smoothed.CNA.object, plot.type="c", cbys.layout=c(2,1), cbys.nchrom=6)
#Plot by plateaus
plot(segment.smoothed.CNA.object, plot.type="p")
```

plotSample

Plot the data and results from segmentation for a single sample

Description

Plots the data for a single sample from a copy number array experiment (aCGH, ROMA etc.) along with the results of segmenting it into regions of equal copy numbers.

Usage

х	an object of class DNAcopy resulting from analyzing data from copy number array experiments.
sampleid	the sample for which the plot is requested. Should be a valid sample name or number. If missing the first sample is plotted.
chromlist	a vector of chromosome numers or names to be plotted. If missing the whole genome is plotted.
xmaploc	a logical indicating if data are plotted against genomic position or Index. De- faults to FALSE.
col	a vector of two colors that can be used for alternating colors for successive chromosomes.
pch	the plotting character. Defaults to
cex	the size of plotting character. If missing it is set to 3 if pch is '.' and 1 otherwise.
altcol	a logical indicating if colors of successive chromosomes should be alternated. Defaults to TRUE.

plotSample

segcol	color for segment means.
zeroline	a logical indicating if the zeroline is drawn. Defaults to TRUE.
zlcol	color for zero line.
lwd	thickness of the lines.
xlab	the x-axis lavel. If missing Index or Genomic Position will be used depending on xmaploc.
ylab	the y-axis label. If missing log(CN) or LOH will be used depending on data type.
main	the main title. If missing sample name will be used.
	other arguments to the plot function can be passed here.

Details

This function plots the whole genome and segmentation results for a single sample. This function overcomes the deficiency in the plot.DNAcopy function which cycles through all the samples. If sampleid is not specified the first sample is plotted.

Examples

#Read in two examples from Snijders et al.

data(coriell)

#Combine into one CNA object to prepare for analysis on Chromosomes 1-23

#We generally recommend smoothing single point outliers before analysis #Make sure to check that the smoothing is proper

smoothed.CNA.object <- smooth.CNA(CNA.object)</pre>

#Segmentation at default parameters

segment.smoothed.CNA.object <- segment(smoothed.CNA.object, verbose=1)</pre>

Plot whole sample c13330

plotSample(segment.smoothed.CNA.object, sampleid="c13330")

Plot only chromosomes 1,3,5,7,9 from first sample
plotSample(segment.smoothed.CNA.object, sampleid=1, chromlist=c(1,3,5,7,9))

segment

Description

This program segments DNA copy number data into regions of estimated equal copy number using circular binary segmentation (CBS).

Usage

Х	an object of class CNA
weights	a vector of weights for the probes. The weights should be inversely proportional to their variances. Currently all weights should be positive i.e. remove probes with zero weight prior to segmentation.
alpha	significance levels for the test to accept change-points.
nperm	number of permutations used for p-value computation.
p.method	method used for p-value computation. For the "perm" method the p-value is based on full permutation. For the "hybrid" method the maximum over the entire region is split into maximum of max over small segments and max over the rest. Approximation is used for the larger segment max. Default is hybrid.
min.width	the minimum number of markers for a changed segment. The default is 2 but can be made larger. Maximum possible value is set at 5 since arbitrary widths can have the undesirable effect of incorrect change-points when a true signal of narrow widths exists.
kmax	the maximum width of smaller segment for permutation in the hybrid method.
nmin	the minimum length of data for which the approximation of maximum statistic is used under the hybrid method. should be larger than 4*kmax
eta	the probability to declare a change conditioned on the permuted statistic exceeding the observed statistic exactly $j (= 1,,nperm*alpha)$ times.
sbdry	the sequential boundary used to stop and declare a change. This boundary is a function of nperm, alpha and eta. It can be obtained using the function "getbdry" and used instead of having the "segment" function compute it every time it is called.
trim	proportion of data to be trimmed for variance calculation for smoothing outliers and undoing splits based on SD.

segment

undo.splits	A character string specifying how change-points are to be undone, if at all. De- fault is "none". Other choices are "prune", which uses a sum of squares criterion, and "sdundo", which undoes splits that are not at least this many SDs apart.
undo.prune	the proportional increase in sum of squares allowed when eliminating splits if undo.splits="prune".
undo.SD	the number of SDs between means to keep a split if undo.splits="sdundo".
verbose	level of verbosity for monitoring the program's progress where 0 produces no printout, 1 prints the current sample, 2 the current chromosome and 3 the current segment. The default level is 1.

Details

This function implements the cicular binary segmentation (CBS) algorithm of Olshen and Venkatraman (2004). Given a set of genomic data, either continuous or binary, the algorithm recursively splits chromosomes into either two or three subsegments based on a maximum t-statistic. A reference distribution, used to decided whether or not to split, is estimated by permutation. Options are given to eliminate splits when the means of adjacent segments are not sufficiently far apart. Note that after the first split the α -levels of the tests for splitting are not unconditional.

We recommend using one of the undoing options to remove change-points detected due to local trends (see the manuscript below for examples of local trends).

Since the segmentation procedure uses a permutation reference distribution, R commands for setting and saving seeds should be used if the user wishes to reproduce the results.

Data that are NA, Inf, NaN will be removed on a per sample basis for "genomdat" and all samples for "chrom" and "maploc".

Value

An object of class DNAcopy. It has three elements:

data	The original CNA object which was the input for segment
out	a data frame with six columns. Each row of the data frame contains a segment for which there are six variables: the sample id, the chromosome number, the map position of the start of the segment, the map position of the end of the segment, the number of markers in the segment, and the average value in the segment.
segRows	a data frame with the start and end row of each segment in the data matrix. print command shows it with the showSegRows=T
call	the call that produced the output object.

Author(s)

Venkatraman E. Seshan <seshanv@mskcc.org> and Adam Olshen <olshena@biostat.ucsf.edu>

References

Olshen, A. B., Venkatraman, E. S., Lucito, R., Wigler, M. (2004). Circular binary segmentation for the analysis of array-based DNA copy number data. *Biostatistics* 5: 557-572.

Venkatraman, E. S., Olshen, A. B. (2007) A faster circular binary segmentation algorithm for the analysis of array CGH data. *Bioinformatics* 23: 657-63.

Examples

```
# test code on an easy data set
set.seed(25)
genomdat <- rnorm(500, sd=0.1) +
rep(c(-0.2,0.1,1,-0.5,0.2,-0.5,0.1,-0.2),c(137,87,17,49,29,52,87,42))
plot(genomdat)
chrom <- rep(1:2,c(290,210))
maploc <- c(1:290,1:210)</pre>
test1 <- segment(CNA(genomdat, chrom, maploc))</pre>
# test code on a noisier and hence more difficult data set
set.seed(51)
genomdat <- rnorm(500, sd=0.2) +</pre>
rep(c(-0.2,0.1,1,-0.5,0.2,-0.5,0.1,-0.2),c(137,87,17,49,29,52,87,42))
plot(genomdat)
chrom <- rep(1:2,c(290,210))
maploc <- c(1:290,1:210)</pre>
test2 <- segment(CNA(genomdat, chrom, maploc))</pre>
# test code for weighted CBS
set.seed(97)
wts <- sample(1:3, 500, replace=TRUE)</pre>
genomdat <- rnorm(500, sd=0.3)/sqrt(wts) +</pre>
rep(c(-0.2,0.1,1,-0.5,0.2,-0.5,0.1,-0.2),c(137,87,17,49,29,52,87,42))
plot(genomdat)
chrom <- rep(1:2,c(290,210))
maploc <- c(1:290,1:210)</pre>
test3 <- segment(CNA(genomdat, chrom, maploc), weights=wts)</pre>
#A real analyis
data(coriell)
#Combine into one CNA object to prepare for analysis on Chromosomes 1-23
CNA.object <- CNA(cbind(coriell$Coriell.05296,coriell$Coriell.13330),</pre>
                   coriell$Chromosome,coriell$Position,
                   data.type="logratio", sampleid=c("c05296", "c13330"))
#We generally recommend smoothing single point outliers before analysis
#Make sure to check that the smoothing is proper
smoothed.CNA.object <- smooth.CNA(CNA.object)</pre>
```

segments.p

segments.p

p-values for the change-points

Description

This program computes pseudo p-values and confidence intervals for the change-points found by the circular binary segmentation (CBS) algorithm.

Usage

```
segments.p(x, ngrid=100, tol=1e-6, alpha=0.05, search.range=100, nperm=1000)
```

х	an object of class DNAcopy
ngrid	number grid points for the integral
tol	tolerance level for calculating nu
alpha	Confidence level is 1-alpha
search.range	statistic is maximized over nu +/- search.range
nperm	number of permutations for confidence interval

Details

The p-values are obtained by applying Siegmund's approximation for the maximal statistic from binary segmenting consecutive segments within a chromosome. This p-value is only to give the relative importance of the change-points as the CBS is different from the algorithm used here.

The confidence intervals are obtained by a permutation algorithm. The data are permuted to the left and right of the identified change-point and the location of the maximal binary segmentation statistic computed. The confidence interval is given by the quantiles of the permutation distribution of the locations.

The statistical properties of this confidence interval is unknown. It is used to give an idea of the uncertainity on the location of the change-points as the CBS is different from the algorithm used here.

Value

a data frame with ten columns. The maximal statistic from binary segmentation, the p-values and lower and upper alpha/2 confidence limits (as genomic positions) are added to the six columns from the segment command.

NOTE: THE p VALUES ARE APPROXIMATE TAIL PROBABILITIES. ANY VALUE GREATER THAN 0.1 CAN HAVE LARGE ERROR. p > 1 ARE REPLACED WITH 1.

Author(s)

Venkatraman E. Seshan

Examples

```
# test code on an easy data set
set.seed(25)
genomdat <- rnorm(500, sd=0.1) +
rep(c(-0.2,0.1,1,-0.5,0.2,-0.5,0.1,-0.2),c(137,87,17,49,29,52,87,42))
plot(genomdat)
chrom <- rep(1:2,c(290,210))
maploc <- c(1:290,1:210)
test1 <- segment(CNA(genomdat, chrom, maploc))
segments.p(test1)
```

segments.summary Additional summary measured for the segments

Description

This program computes the standard deviation, median and the mad of the data for each segment found by the CBS algorithm.

Usage

segments.summary(x)

16

smooth.CNA

Arguments

х

an object of class DNAcopy

Value

a data frame with nine columns. The sd, median and mad of each segment is added to the six columns from the segment command.

Author(s)

Venkatraman E. Seshan

Examples

```
# test code on an easy data set
set.seed(25)
genomdat1 <- rnorm(500, sd=0.1) +
rep(c(-0.2,0.1,1,-0.5,0.2,-0.5,0.1,-0.2),c(137,87,17,49,29,52,87,42))
genomdat2 <- rnorm(500, sd=0.1) +
rep(c(-0.2,0.1,1,-0.5,0.2,-0.5,0.1,-0.2),c(137,87,17,49,29,52,87,42))
genomdat1[sample(1:500,5)] <- NA
chrom <- rep(1:2,c(290,210))
maploc <- c(1:290,1:210)
test1 <- segment(CNA(cbind(genomdat1,genomdat2), chrom, maploc))
segments.summary(test1)
```

smooth.CNA	Smooth a 'Copy Number Array' data object	

Description

Detect outliers and smooth the data prior to analysis by programs such as circular binary segmentation (CBS).

Usage

х	Copy number array data object
smooth.region	number of points to consider on the left and the right of a point to detect it as an outlier. (default=10)
outlier.SD.scale	
	the number of SDs away from the nearest point in the smoothing region to call a point an outlier.

smooth.SD.scale	
	the number of SDs from the median in the smoothing region where a smoothed point is positioned.
	proportion of data to be trimmed for variance calculation for smoothing outliers and undoing splits based on SD.

Value

An object of class CNA with outliers smoothed i.e the logratio values of singleton outliers is shrunk towards the values of its neighbors. The output is of the same dimension as the input.

Examples

data(coriell)

#Combine into one CNA object to prepare for analysis on Chromosomes 1-23

#We generally recommend smoothing single point outliers before analysis #Make sure to check that the smoothing is proper

```
smoothed.CNA.object <- smooth.CNA(CNA.object)</pre>
```

subset.CNA

Subset a 'Copy Number Array' data object

Description

Function to return a subset of a copy number array data object by a list of chromosomes and sample.

Usage

```
## S3 method for class 'CNA'
subset(x, chromlist=NULL, samplelist=NULL, ...)
```

х	Copy number array data object
chromlist	chromosomes of interest. Should be a subset of the valid chromosome names in the original data.
samplelist	samples of interest. Can be integers denoting the samples of interest or a vector of valid sample names.
	other arguments which may be passed to subset.

subset.DNAcopy

Value

An object of class CNA with the data for the list of chromosomes and samples of interest.

Examples

```
#subset.CNA.object <- subset.CNA(CNA.object,chromlist=1:10,samplelist="c05296")
subset.CNA.object <- subset(CNA.object,chromlist=1:10,samplelist="c05296")</pre>
```

subset.DNAcopy

Subset a DNAcopy data object

Description

Function to return a subset of a copy number array data object by a list of chromosomes and sample.

Usage

S3 method for class 'DNAcopy'
subset(x, chromlist=NULL, samplelist=NULL, ...)

Arguments

х	DNAcopy object
chromlist	chromosomes of interest. Should be a subset of the valid chromosome names in the original data.
samplelist	samples of interest. Can be integers denoting the samples of interest or a vector of valid sample names.
	other arguments which may be passed to subset.

Value

An object of class DNAcopy with the input data and the results of segmenting them only for the chromosomes and samples of interest.

zoomIntoRegion

Description

This program computes the frequency of gains and losses for each probe as a function of level of mad.

Usage

```
zoomIntoRegion(x, chrom, sampleid, maploc.start=NULL, maploc.end=NULL,
    pt.pch=NULL, pt.cex=NULL, pt.col=NULL, segcol=NULL, seglwd=NULL,
    main=NULL, xlab=NULL, ylab=NULL, ...)
```

Arguments

х	an object of class DNAcopy.
chrom	the chromosome in which the region lies.
sampleid	the sample of interest.
maploc.start	genomic start position of the region of interest. Default is the beginning of the chromosome.
maploc.end	genomic end position of the region of interest. Default is the end of the chromo- some.
pt.pch	the plotting character used for plotting the log-ratio values (default is ".").
pt.cex	the size of plotting character used for the log-ratio values (default is 3 if "." and 1 otherwise).
pt.col	the color used for the points. Default is green3.
segcol	the color of the lines indicating the segment means. If missing the line color is set to be red.
seglwd	line weight of lines for segment mean and zeroline. If missing it is set to 3.
main	figure title. If missing will be generated by pasting the chromosome, range and sample name together.
xlab	x-axis label. If missing "Genomic position" will be used
ylab	y-axis label. If missing "log-ratio" will be used
•••	additional plotting options.

Details

This command plots the region of interest with the log-ratio and segments. It works for a region from a single chromosome in a single sample. So if more than one chromosome and/or one sample are given only the first chromosome from the first sample will be used.

zoomIntoRegion

Author(s)

Venkatraman E. Seshan <seshanv@mskcc.org>

Examples

data(coriell)

#Combine into one CNA object to prepare for analysis on Chromosomes 1-23

#We generally recommend smoothing single point outliers before analysis #Make sure to check that the smoothing is proper

smoothed.CNA.object <- smooth.CNA(CNA.object)</pre>

#Segmentation at default parameters

segment.smoothed.CNA.object <- segment(smoothed.CNA.object, verbose=1)</pre>

zoomIntoRegion(segment.smoothed.CNA.object, chrom=10, sampleid="c05296")

Index

* datasets coriell, 3 cytoBand, 4 * internal DNAcopy-internal, 5 * nonparametric CNA, 2DNAcopy, 4 exon.segment, 6 getbdry, 7 glFrequency, 7 plot.DNAcopy, 8 plotSample, 10 segment, 12segments.p, 15 segments.summary, 16 smooth.CNA, 17 subset.CNA, 18 subset.DNAcopy, 19 zoomIntoRegion, 20 changepoints (DNAcopy-internal), 5 CNA, 2coriell, 3 cytoBand, 4 default.DNAcopy.bdry(getbdry),7 DNAcopy, 4 DNAcopy-internal, 5 exon.changepoint (DNAcopy-internal), 5 exon.segment, 6 getbdry, 7 glFrequency, 7inflfact (DNAcopy-internal), 5 plot.DNAcopy, 8 plotSample, 10 print.CNA(CNA), 2

print.DNAcopy(DNAcopy), 4

segment, 12
segments.p, 15
segments.summary, 16
smooth.CNA, 17
smooth.data(DNAcopy-internal), 5
subset.CNA, 18
subset.DNAcopy, 19

trimmed.variance (DNAcopy-internal), 5

zoomIntoRegion, 20